



In London, 7.5.2019

## LAURA: Calibration procedure

To Whom It May Concern:

We Erba Lachema s.r.o. with registered office in Czech Republic, Karásek 2219/1D, 621 00 Brno, represented by its Product Manager Veronika Plešková hereby confirm as producer of LAURA instrument that this device is firmly calibrated by manufacturer and does not require any regular calibration by the user. Correct calibration of the instrument is controlled by grey strips provided initially together with the instrument.

  
Veronika Plešková

  
Erba Corporate Services Ltd.  
GW1 Great West House  
Great West Road, London TW8 9DF

Product manager, Urine Analysis

Erba Corporate Service Ltd.  
GW1 Great West House  
Great West Road, London TW8 9DF



DEKAPHAN LAURA

Seq.No: 0001

URINORM N

2022.08.13

11:39 AM

.....  
BLD NEG  
BIL NEG  
UBG NORM  
KET NEG  
GLU NEG  
PRO NEG  
PH 6  
NIT NEG  
LEU NEG  
SG 1.020

*Ahmad*  
13/8/22

QC TEST: OK

DEKAPHAN LAURA

Seq.No: 0002

URINORM P

2022.08.13

11:40 AM

.....  
\*BLD 50 Ery/u $\bar{l}$   
\*BIL 6 mg/d $\bar{l}$   
\*UBG 3 mg/d $\bar{l}$   
\*KET 156 mg/d $\bar{l}$   
\*GLU 500 mg/d $\bar{l}$   
\*PRO 100 mg/d $\bar{l}$   
PH 9  
\*NIT POS  
\*LEU 500 Leu/u $\bar{l}$   
SG 1.000

*Ahmad*  
13/8/22

QC TEST: OK