




In London, 7.5.2019

LAURA: Calibration procedure

To Whom It May Concern:

We Erba Lachema s.r.o. with registered office in Czech Republic, Karásek 2219/1D, 621 00 Brno, represented by its Product Manager Veronika Plešková hereby confirm as producer of LAURA instrument that this device is firmly calibrated by manufacturer and does not require any regular calibration by the user. Correct calibration of the instrument is controlled by grey strips provided initially together with the instrument.

Veronika Plešková



Product manager, Urine Analysis

Erba Corporate Service Ltd.
GW1 Great West House
Great West Road, London TW8 9DF



14/8/22

2022.08.14 09:02

DEKAPHAN LAURA
Seq.No: 0001
URINORM N
2022.08.14 09:53

BLD NEG
BIL NEG
UBG NORM
KET NEG
GLU NEG
PRO NEG
PH 6
NIT NEG
LEU NEG
SG 1.030

QC TEST: OK

DEKAPHAN LAURA
Seq.No: 0002
URINORM P
2022.08.14 09:56

*BLD 250 Ery/u1
*BIL 6 mg/d1
*UBG 3 mg/d1
*KET 52 mg/d1
*GLU 250 mg/d1
*PRO 500 mg/d1
pH 9
*NIT POS
*LEU 125 Leu/u1
SG 1.000

QC TEST: OK