



March 09, 2022

*Calibration Certificate for Automated Clinical Chemistry Analyzer
Model: XL 640 (Sl. No. 640371)*

This is to certify that the Automated Clinical Chemistry Analyzer, Model: XL 640 installed at Redcliffe Lifetech Pvt. Ltd., 2nd Floor, M.P. Jewellery, Lalganesh, Guwahati, Assam, PIN - 781034 has been calibrated. All parameters have been checked and found well within the limit. Machine precision and all parameters control found within range.

This calibration is valid till March 08, 2023.

For Transasia Bio-Medicals Ltd.




*Taraknath Chakraborty,
Dy. Regional Service Manager,
Guwahati, Assam*

Result Reprint

Report Type : Controls 17-Sep-2022

Consumable: RANDOX L2

Gr #	Lot #	Consumable	Test	Result Unit	Flag	Curve #	Result Date	Mean	SD	Interval (3SD)
1	1543UN	RANDOX L2	Na	135.50 mmol/l		0	17-Sep-2022 11:29:10	138.000	3.500	127.5 - 148.5
2	1543UN	RANDOX L2	K	3.65 mmol/l	-1SD	0	17-Sep-2022 11:29:10	3.900	0.150	3.45 - 4.35
3	1543UN	RANDOX L2	Cl	92.30 mmol/l		0	17-Sep-2022 11:29:10	96.000	3.920	84.24 - 107.76
4	1543UN	RANDOX L2	SGPTD	39.7 U/L		146012	17-Sep-2022 11:31:55	41.000	4.000	29 - 53
5	1543UN	RANDOX L2	SGOTD	33.3 U/L		146014	17-Sep-2022 11:32:13	36.000	3.500	25.5 - 46.5
6	1543UN	RANDOX L2	GGT	54.9 U/L		146018	17-Sep-2022 11:32:49	53.223	4.000	41.223 - 65.223
7	1543UN	RANDOX L2	AMYY	69 U/L		146020	17-Sep-2022 11:33:07	66.000	5.000	51 - 81
8	1543UN	RANDOX L2	LIP	47.9 U/L	+1SD	146022	17-Sep-2022 11:33:25	40.000	4.000	28 - 52
9	1543UN	RANDOX L2	PHOS	4.34 mg/dl		146025	17-Sep-2022 11:33:52	4.060	0.310	3.13 - 4.99
10	1543UN	RANDOX L2	BIDD	0.91 mg/dl		146027	17-Sep-2022 11:34:10	0.916	0.140	0.496 - 1.336
11	1543UN	RANDOX L2	BIT	1.64 mg/dl		146029	17-Sep-2022 11:34:28	1.680	0.170	1.17 - 2.19
12	1543UN	RANDOX L2	UREA	37.7 mg/dl	-1SD	146031	17-Sep-2022 11:34:46	43.900	3.300	34 - 53.8
13	1543UN	RANDOX L2	CRENZ	1.38 mg/dl		146033	17-Sep-2022 11:35:04	1.470	0.140	1.05 - 1.89
14	1543UN	RANDOX L2	GLU	119.2 mg/dl		146035	17-Sep-2022 11:35:23	111.000	8.400	85.8 - 136.2
15	1543UN	RANDOX L2	TRIG	95.5 mg/dl		146037	17-Sep-2022 11:35:41	96.500	7.650	73.55 - 119.45
16	1543UN	RANDOX L2	CHOL	165 mg/dl		146039	17-Sep-2022 11:35:59	166.467	10.500	134.967 - 197.967
17	1543UN	RANDOX L2	HDLc	49.3 mg/dl		146042	17-Sep-2022 11:36:26	52.500	3.850	40.95 - 64.05

Print Date 17-Sep-2022

Result Reprint

Report Type : Controls 17-Sep-2022

Consumable: RANDOX L2

Run	Reagent	Unit	Result	Curve #	Run Date	Mean	SD	Interval (3SD)
18	RANDOX L2	UA	5.6 mg/dl	146044	17-Sep-2022 11:36:44	5.750	0.370	4.64 - 6.86
19	RANDOX L2	CA	9.5 mg/dl +1SD	146046	17-Sep-2022 11:37:02	9.079	0.420	7.819 - 10.339
20	RANDOX L2	ALDO	3.72 g/dl	146048	17-Sep-2022 11:37:20	3.613	0.310	2.883 - 4.743
21	RANDOX L2	PRO	5.10 g/dl -1SD	146050	17-Sep-2022 11:37:38	5.740	0.570	4.03 - 7.45
22	RANDOX L2	FEE	121.4 µg/dl +1SD.W*	146061	17-Sep-2022 12:11:53	103.000	9.300	75.1 - 130.9
23	RANDOX L2	CA	8.3 mg/dl -1SD	146081	17-Sep-2022 12:58:15	9.079	0.420	7.819 - 10.339
24	RANDOX L2	ALPU	212 U/L	146299	17-Sep-2022 14:40:56	203.000	15.000	158 - 248

Result Reprint

Report Type : Controls 17-Sep-2022 Consumable: RANDOX L3

Sl#	Lot #	Consumable	Test	Result Unit	Flag	Curve #	Result Date	Mean	SD	Interval (3SD)
1	1211UE	RANDOX L3	Na	148.30 mmol/l	-1SD	0	17-Sep-2022 11:29:57	154.000	3.500	143.5 - 164.5
2	1211UE	RANDOX L3	K	5.56 mmol/l	-1SD	0	17-Sep-2022 11:29:57	5.950	0.240	5.23 - 6.67
3	1211UE	RANDOX L3	Cl	110.60 mmol/l		0	17-Sep-2022 11:29:57	111.000	4.500	97.5 - 124.5
4	1211UE	RANDOX L3	SGPTD	112.9 U/L	-1SD	146013	17-Sep-2022 11:32:04	127.000	12.500	89.5 - 164.5
5	1211UE	RANDOX L3	SGOTD	139.3 U/L		146015	17-Sep-2022 11:32:22	136.000	13.500	95.5 - 176.5
6	1211UE	RANDOX L3	GGT	194.9 U/L		146019	17-Sep-2022 11:32:58	187.345	13.000	148.345 - 226.345
7	1211UE	RANDOX L3	AMYY	260 U/L		146021	17-Sep-2022 11:33:16	254.000	19.000	197 - 311
8	1211UE	RANDOX L3	LIP	58.7 U/L		146023	17-Sep-2022 11:33:34	58.000	5.500	41.5 - 74.5
9	1211UE	RANDOX L3	PHOS	7.22 mg/dl		146026	17-Sep-2022 11:34:01	7.220	0.540	5.6 - 8.84
10	1211UE	RANDOX L3	BIDD	1.94 mg/dl		146028	17-Sep-2022 11:34:19	1.970	0.210	1.34 - 2.6
11	1211UE	RANDOX L3	BIT	4.98 mg/dl		146030	17-Sep-2022 11:34:37	5.180	0.550	3.53 - 6.83
12	1211UE	RANDOX L3	UREA	106.3 mg/dl	-1SD	146032	17-Sep-2022 11:34:55	118.000	9.000	91 - 145
13	1211UE	RANDOX L3	CRENZ	4.18 mg/dl		146034	17-Sep-2022 11:35:13	4.240	0.430	2.95 - 5.53
14	1211UE	RANDOX L3	GLU	304.1 mg/dl	+1SD	146036	17-Sep-2022 11:35:31	276.000	21.000	213 - 339
15	1211UE	RANDOX L3	TRIG	259.5 mg/dl		146038	17-Sep-2022 11:35:50	267.000	21.000	204 - 330
16	1211UE	RANDOX L3	CHOL	293 mg/dl		146040	17-Sep-2022 11:36:08	281.000	18.000	227 - 335
17	1211UE	RANDOX L3	HDLc	88.0 mg/dl		146043	17-Sep-2022 11:36:35	85.300	9.000	58.3 - 112.3

Print Date 17-Sep-2022

Result Reprint

Report Type : Controls 17-Sep-2022

Consumable: RANDOX L3

Sr#	Lot #	Consumable	Test	Result Unit	Flag	Curve #	Result Date	Mean	SD	Interval (3SD)
18	1211UE	RANDOX L3	UA	8.9 mg/dl		146045	17-Sep-2022 11:36:53	9.260	0.610	7.43 - 11.09
19	1211UE	RANDOX L3	ALBD	3.00 g/dl		146049	17-Sep-2022 11:37:29	3.050	0.230	2.36 - 3.74
20	1211UE	RANDOX L3	PRO	4.59 g/dl		146051	17-Sep-2022 11:37:47	4.590	0.460	3.21 - 5.97
21	1211UE	RANDOX L3	FEE	230.4 ug/dl		146062	17-Sep-2022 12:12:02	214.837	19.000	157.837 - 271.837
22	1211UE	RANDOX L3	CA	12.7 mg/dl		146085	17-Sep-2022 13:18:17	12.700	0.650	10.75 - 14.65
23	1211UE	RANDOX L3	ALPU	351 U/L	-1SD	146330	17-Sep-2022 15:01:09	386.000	29.000	299 - 473

Result Reprint

Report Type : Controls 17-Sep-2022

Consumable: HbA1C CONTROL H

Sr #	Lot #	Consumable	Test	Result Unit	Flag	Curve #	Result Date	Mean	SD	Interval (3SD)
1	S012245	HbA1C CONTROL H	HA1cD	14.3 %		146058	17-Sep-2022 12:11:26	13.400	1.000	10.4 - 16.4

Result Reprint

Report Type : Controls

17-Sep-2022

Consumable: HBA1C LOW

Sr #	Lot #	Consumable	Test	Result Unit	Flag	Curve #	Result Date	Mean	SD	Interval (3SD)
1	S022219	HBA1C LOW	HA1GD	5.8 %		146056	17-Sep-2022 11:54:10	5.500	0.400	4.3 - 6.7

