RML EQAS CORRECTIVE & PREVENTIVE ACTION

Name of the department: Clinical Biochemistry

EQAS Programme

: RML EQAS

Cycle / Sample No

: 11

Month and year

: AUGUST - 2022

Nature of Non conformance:

Serum SGOT and Bilirubin SDI: ≥ ±2.0

Corrective and preventive action taken:

Bilirubin and SGOT Reagent stability, Expiry status and internal quality Control values were verified on the day (12/07/22) we processed our EQAS sample and the QC value is within the normal range only. It may be due to Random Error.

Review in the next EQAS sample.

Date:

Meln

Approved by:

Dr Mohd Khaleel Msc.Ph.D Head-Clinical Laboratory