

Brno, 10.11.2016

**Laura Smart: Calibration procedure**

To Whom It May Concern:

We Erba Lachema s.r.o. with registered office in Czech Republic, Karásek 1d, 621 00 Brno, represented by its Global Product Manager Pavel Suchan hereby confirm as producer of the Laura Smart instrument that this device is firmly calibrated by manufacturer and does not require any regular calibration by the user. Correct calibration of the instrument is controlled by grey strips provided initially together with the instrument.

*P. Suchan*  
RNDr. Pavel Suchan Ph.D.

**Erba Lachema s.r.o.**  
Karásek 1d, 621 00 Brno  
DIČ: CZ26918846 -10-

Global Product Manager, Urine Analysis

Erba Lachema s.r.o.

Karásek 1d, Brno 62100, Czech Republic



DEKAPHAN LAURA  
Seq.No: 0001  
URINORM N  
2022.11.04

09:20

BLD NEG  
BIL NEG  
UBG NORM  
KET NEG  
GLU NEG  
PRO NEG  
pH 6  
NIT NEG  
LEU NEG  
SG 1.025

QC TEST: OK  
<Admin>

DEKAPHAN LAURA  
Seq.No: 0002  
URINORM P  
2022.11.04

09:21

\*BLD 10 Ery/u1 1+  
\*BIL 3 mg/dl 2+  
\*UBG 1 mg/dl 1+  
\*KET 52 mg/dl 2+  
\*GLU 1000 mg/dl 3+  
\*PRO 100 mg/dl 2+  
pH 9  
\*NIT POS  
\*LEU 75 Leu/u1 1+  
SG 1.000

QC TEST: OK  
<Admin>

*Abhishek*  
*3/11/22*