

Date: 21.12.2021

To, Right Care Pathology Lab Ms. Jyoti D Phadtare H. No.3823/1. S. No 31 Ganesh Nagar, Dhankawadi Pune Maharashtra 411 043.

Dear Sir,

This is to certify that, the Instrument Model AGD300 Fully Automatic Hemaotolgy Analyzer bearing Serial No. 0740619193650 has been Calibrated and checked with CBC- 3D Normal control print out is given.

Kindly find attached service report of our Engineer for your reference.

Thanking you and assuring you of our attention & prompt services at all times.

Sincerely,

For AGD ALLEGICALS (P) Ltd.

Authorized Signatory



### AGD BIOMEDICALS PVT. LTD.

# CBC-3D

#### HEMATOLOGY CONTROL

LOT: B0822 EXP: 2022/11/05

1	Low		mal	Hi	High		
	Parameter	Assay Mean	Range (+/-)	Assay Mean	Range (+/-)	Assay Mean	Range (+/-)
	WBC	2.0	± 0.5	8.1	± 1.0	20.9	± 2.5
	LYM (%)	62.0	± 12.0	30.9	± 8.0	12.5	± 6.0
	MID (%)	5.8	± 5.8	5.3	± 5.3	3.8	± 3.8
	NEUT (%)	32.2	± 10.0	63.8	± 8.0	83.7	± 8.0
	LYM (#)	1.3	± 0.3	2.5	± 0.7	2.6	±1.3
	MID (#)	0.1	± 0.1	0.4	± 0.4	0.8	± 0.8
AGD	NEUT (#)	0.7	± 0.3	5.1	± 0.7	17.4	± 1.7
300	RBC	2.13	± 0.18	4.22	± 0.24	5.17	± 0.30
500	HGB	59	± 4	139	±6	191	± 8
	нст	16.8	± 1.5	37.9	± 2.0	50.9	± 2.4
	MCV	79.1	± 5.0	89.9	± 5.0	98.4	± 5.0
	мсн	27.4	± 2.5	32.8	± 2.5	36.8	± 2.5
	MCHC	348	± 30	365	± 30	375	± 30
	RDW-CV	15.7	± 3.0	14.7	± 3.0	14.3	± 3.0
	PLT	72	± 20	259	± 40	498	± 60
	MPV	9.8	± 3.0	10.1	± 3.0	9.8	±.3.0

Note: Assay value provided by AGD Biomedicals Pst. Ltd.

Exclusive Supplier of R&D SYSTEMS Hematology controls



**Installation Qualification** 

Instrument Type: Blood Cell Counter	Manufacture Name: AGD Biomedicals (P) LTD
Model Name: AGD 300	Supplier Name: AGD Biomedicals (P) LTD

**Operational Verification Report** 

Instrument Type: Blood Cell Counter	Name of the Customer: Right Care Pathology Lab	
Serial Number: 0740619193650	Model Number: AGD 300	
Date of Installation: 05.08.2019	Date of receipt: 04.08.2019	

**Power Check:** 

Power Supply:	220 V AC	
Earth Voltage	OK	

Physical Check of the equipment

Check for damage	Checked	No damage found
Visual & functional inspection	Checked	OK
No damage on internal assembly	Checked	No damage found

System configuration equipment

Found OK
and the second
Found OK
No Interference

#### **Result:**

Installed the AGD 300 satisfactorily and is ready for operational verification

Installed by: Mr Mukund T	Checked by: Bhyshan Grampurit
Signature:	Signature: Slundy
Date: 05/08/	2019. Date: 05/08/2019
7	R   9   .

E: sales@agdbio.com
W: www.agdbio.com



#### **Operational Qualification**

Instrument Type: Fully Automatic	Manufacture Name: AGD Biomedicals
Hemaotolgy Analyzer	(P) LTD
Model Name: AGD 300	Supplier Name: AGD Biomedicals (P) LTD
Installed On: 05.08.2019	Service Engineer: Mr Mukund T

#### **Operational Verification Report**

Name of the Customer: Right Care Pathology Lab	Instrument Type: Fully Automatic Hemaotolgy Analyzer	
Model Number: AGD 300	Serial Number: 07460919193650	
<b>Date of receipt: 04.08.2019</b>	Date of Installation: 05.08.2019	

### **Operational Check:**

Background Check	OK	
Control Check	All Values are within Range	
Reproducibility Check	OK	

Result:

Operational verification of AGD 300 has been completed satisfactorily and it is ready for the performance verification

Operations performed by: Mr. Mukund

T

Signature:

Date:

05/08/2019.

Checked by:

Bhushan Grampwohit

Signature:

Date: 05/08/2019.



## Performance Qualification

Instrument Type: Fully Automatic Hemaotolgy Analyzer	Manufacture Name: AGD Biomedicals (P) LTD
Model Name: AGD 300	Supplier Name: AGD Biomedicals (P) LTD
Installed On: 05.08.2019	Service Engineer: Mr Mukund T

**Performance Verification Report** 

Name of the Customer: Right Care Pathology Lab	Instrument Type: Fully Automatic Hemaotolgy Analyzer
Model Number: AGD 300	Serial Number: 07460919193650
Date of receipt: 04.08.2019	Date of Installation: 05.08.2019

**Performance Check:** 

Performance Check.	
Background Check	OK
Pressure Check	OK
HGB Unit Check	OK
Syringe Movement Check	OK
Printer Check	OK

E : sales@agdbio.com

W: www.agdbio.com



Result:

The performance verification of AGD 300 has been completed satisfactorily and

Instrument is ready to be placed in service.

	HOD: Bhushen Grangawiit
) A C ( ) ( ) ( ) ( )	Signature: Sheety Date: 05/08/2019
Mr. Sohan Rathod.	Date: 05/08/2019
Performance check:	Checked by:
Performed by: Mr. Mukund Tmape	Mr. Mukund Jonuse
Date: 5108/2019	Date: 510812019.

Installation, Operational and Performance Qualification Certificate

The equipment AGD 300 Serial No 0740619193650 is installed and commissioned. The operational and performance verification is done and observations are:

- All system components are working in normal condition
- All tests were carried out in accordance with qualification protocol and standard operating procedure.
- The performance of the AGD 300 is as per the acceptance criteria.



Date: 23.08.2022

To,
Right Care Pathology Lab
Ms. Jyoti D Phadtare
I. No.3823/1. S. No 31
Ganesh Nagar, Dhankawadi
Pune
Maharashtra 411 043.

Dear Sir,

This is to certify that, the Instrument Model AGD300 Fully Automatic Hemaotolgy Analyzer bearing Serial No. 0740619193650 has been Calibrated and checked with CBC- 3D Normal control print out is given.

Kindly find attached service report of our Engineer for your reference.

Thanking you and assuring you of our attention & prompt services at all times.

Sincerely,

For AGD Biomedicals (P) Ltd.

**Authorized Signatory**