



In London, 7.5.2019

## LAURA: Calibration procedure

To Whom It May Concern:

We Erba Lachema s.r.o. with registered office in Czech Republic, Karásek 2219/1D, 621 00 Brno, represented by its Product Manager Veronika Plešková hereby confirm as producer of LAURA instrument that this device is firmly calibrated by manufacturer and does not require any regular calibration by the user. Correct calibration of the instrument is controlled by grey strips provided initially together with the instrument.

  
Veronika Plešková

  
Erba Corporate Services Ltd.  
GW1 Great West House  
Great West Road, London TW8 9DF

Product manager, Urine Analysis

Erba Corporate Service Ltd.  
GW1 Great West House  
Great West Road, London TW8 9DF





LAURA  
\*\*\*\*\*  
28.12.2022 11:02 AM

DEKAPHAN LAURA  
Seq.No: 0001  
URINORM N  
28.12.2022 11:03 AM

BLD NEG  
BIL NEG  
UBG NORM  
KET NEG  
GLU NEG  
PRO NEG  
PH 6  
NIT NEG  
LEU NEG  
SG 1.020

QC TEST: OK

DEKAPHAN LAURA  
Seq.No: 0002  
URINORM P  
28.12.2022 11:04 AM

\*BLD 10 Ery/u1  
\*BIL 3 mg/dl  
\*UBG 1 mg/dl  
\*KET 52 mg/dl  
\*GLU 500 mg/dl  
\*PRO 100 mg/dl  
pH 9  
\*NIT POS  
\*LEU 75 Leu/u1  
SG 1.000

QC TEST: OK

*Rec by Kusuma  
28/12/22  
Sudhar  
28/12/22*