

**CERTIFICATE OF CALIBRATION**

ADL/MO/SER/018/23-24

**Instrument : Mispa Nano Fully Automated Biochemistry Analyser****Customer Name : Disha CRL Laboratory, Parel****Date of installation : 14/01/2018****Serial No : WN-71106638****Alignment Check**

Motion Check	Checked	Remark
Sample Probe	<input checked="" type="checkbox"/>	Ok
Reagent Probe	<input checked="" type="checkbox"/>	Ok
Mixing Bar	<input checked="" type="checkbox"/>	Ok
Reaction Disk	<input checked="" type="checkbox"/>	Ok
Reagent Disk	<input checked="" type="checkbox"/>	Ok
Sample Disk	<input checked="" type="checkbox"/>	Ok
Fluidic disk	<input checked="" type="checkbox"/>	Ok
Lamp	<input checked="" type="checkbox"/>	Ok

System	Checked	Remark
Main unit	<input checked="" type="checkbox"/>	Ok
Temp. unit	<input checked="" type="checkbox"/>	Ok
Mixing unit	<input checked="" type="checkbox"/>	Ok
Reaction unit	<input checked="" type="checkbox"/>	Ok
Reagent unit	<input checked="" type="checkbox"/>	Ok
Sample unit	<input checked="" type="checkbox"/>	Ok

**AGAPPE DIAGNOSTICS LTD.** ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC15413CORPORATE OFFICE / REAGENT PLANT  
Agappe Hills, Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India.  
Tel: +91 484 286 7000 | Email: agappe@agappe.inEQUIPMENT PLANT  
X/588-CB, Block No. 32, KINFRA Small Industrial Park,  
Melad, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.KOLKATA OFFICE  
406, Merin Matrix, Plot No-10, Block-DN, Sector V,  
Salt Lake City, Kolkata - 700 091.  
Tel: +91 33 4003 0451 | Email: kolkatoffice@agappe.inMUMBAI [REGISTERED OFFICE]  
401 & 402 4<sup>th</sup> Floor, Jaicinch Business Centre, 119DELHI OFFICE  
DSM 540 5<sup>th</sup> Floor, DLF Tower, Shivaji MargBANGALORE OFFICE  
5-6, 11<sup>th</sup> Floor, Red Cross Bhavan, No:26 Race Course Road.



**OPTICAL SECTION**

Filter	Dark Current	Background
340	148	27 339
405	120	31 356
450	107	34 348
510	105	37 803
546	104	34 851
578	104	38 024
630	102	39 948
670	99	38 474
Reference		

Power Supply : 5V  12V  24V

The equipment is calibrated for the different parameters by using Multi Calibrator (Randox CAL 3)  
 Lot.No: 1262UE Exp: 28/01/2025.

This is to certify that the above mentioned product has been calibrated according to the standard operating procedure given by the MINDRAY Medical Systems.

Date of Calibration : **05-04-2023**  
 Next Due : **03-10-2023**

**For Agappe Diagnostics Ltd.**

**Authorised signatory**



**Name: Sahajanand Govind Panchal**  
**Designation: Senior Engineer – Customer Service**

**AGAPPE DIAGNOSTICS LTD.** ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN : U24239MH1998PLC115413

**CORPORATE OFFICE / REAGENT PLANT**  
 Agappe Hills, Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India.  
 Tel. +91 484 286 7000 | Email: agappe@agappe.in

**EQUIPMENT PLANT**  
 X/588 CB, Block No. 32, KINFRA Small Industrial Park,  
 Neilad, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

**KOLKATA OFFICE**  
 406, Merlin Matrix, Plot No-10, Block-DN, Sector V,  
 Salt Lake City, Kolkata - 700 091.  
 Tel: +91 33 4003 0451 | Email: kolkataoffice@agappe.in

**MUMBAI (REGISTERED OFFICE)**  
 401 & 402, 4<sup>th</sup> Floor, Jaisingh Business Centre, 119,

**DELHI OFFICE**  
 DSM 540, 5<sup>th</sup> Floor, DLF Tower, Shivaji Marg,

**BANGALORE OFFICE**  
 S-6, II<sup>nd</sup> Floor, Red Cross Bhavan, No.26 Race Course Road,