

Ortho Clinical Diagnostics

403, Leela Business Park, Andheri Kurla Road,
Andheri East, Mumbai - 400059
T : +91 22 6787 9300
F : +91 22 6787 9333

Calibration Certificate

The below mentioned instrument has been calibrated and tests performed to check the system performance.

Instrument : VITROS XT7600
Serial No : 7600 0983
Customer Name : USHAKAL ABHINAV SPECIALITY HOSPITAL, SANGLI
Calibration performed on : 28/01/2023

The system's calibration includes optics calibration and checking the reproducibility performance of the instrument as per the guidelines provided by the manufacturer.

Next Calibration will be performed in **JULY 2023**.

For **Ortho Clinical Diagnostics India Pvt Ltd.**



Onkar Gajanan Deshpande
Field Engineer - Ortho Care Service
Pune

Date: 28/01/2023

Diagnostics Mode
Assays completed

22/2/2023

11: 21: 16

V3. 7. 2

620301



Status



Samples



Results



QC



Reagents



Diagnostics



Options



Docs



Conditions



1



DIAGNOSTICS - REFLECTOMETER Adjustments

Read Sync By Slot Position (msec/2)

Slot	Tangential	Radial	Slot	Tangential	Radial	Slot	Tangential	Radial
1.	400	256	13.	398	255	25.	400	256
2.	402	255	14.	397	256	26.	401	255
3.	401	255	15.	399	256	27.	402	255
4.	401	255	16.	396	256	28.	402	255
5.	405	255	17.	398	256	29.	402	255
6.	405	255	18.	393	256	30.	403	255
7.	403	255	19.	398	255	31.	403	256
8.	401	255	20.	398	255	32.	399	256
9.	401	255	21.	398	255	33.	401	255
10.	398	255	22.	399	255	34.	401	256
11.	399	255	23.	400	256	35.	401	255
12.	398	255	24.	401	255	36.	388	253

Auto-adjustment complete; touch 'Save' to use the new value or 'Stop' to discard it.



Return



Save



Start



Stop



Correction Factors



Help

Diagnostics Mode
Assays completed

22/2/2023

11:30:24

V3.7.2

620301



Status



Samples



Results



QC



Reagents



Diagnostics



Options



Docs



Conditions



3



A



Keyboard

DIAGNOSTICS - Reflectometer Performance Tests

Tests:

Dark CCD & PD

Static

Static Interleaved

Reference Dynamic

Ref Dynamic XT

MTF

Slot

1

Filter	Cutoff Value
340	17.0
400	14.6
460	15.6
540	16.4
600	16.6
630	16.6
670	17.7
680	17.9

Performance test completed successfully.



Return



Start



Stop



Help

Diagnostics Mode
Assays completed

28/1/2023

18:27:21

V3.7.2

620301



Status



Samples



Results



QC



Reagents



Diagnostics



Options



Docs



Conditions



DIAGNOSTICS - Auto Centering Adjustment

Slide Type:

ST Slides

XT Slides

Location	Dispense Blade		Metering Truck	
	Current	New	Current	New
XT Spot 1	28864	28856	-459675	-459485
ST Spot	28048	28048	-452875	-452723
XT Spot 2	27232	27240	-446075	-445961

Auto-adjustment successful. No Update required. Select 'Return' button to retain current values or 'Save' to update.



Return



Save



Start



Stop



Help

Diagnostics Mode
Assays completed

Navigation icons: Status, Samples, Results, QC, Reagents, Diagnostics, Options, V-Docs, Conditions, and a numeric keypad.

DIAGNOSTICS - Auto Centering Adjustment

Slide Type: ST Slides XT Slides

Location	Dispense Blade		Metering Truck	
	Current	New	Current	New
XT Spot 1	28856	28800	-459381	-459343
ST Spot	28064	28028	-452648	-452705
XT Spot 2	27272	27256	-445915	-446067

Auto-adjustment successful. No Update required. Select 'Return' button to retain current values or 'Save' to update.

Control buttons: Return, Save, Start, Stop, Help