



CHRISTIAN MEDICAL COLLEGE

DEPARTMENT OF CLINICAL BIOCHEMISTRY
CMC EXTERNAL QUALITY ASSURANCE SCHEME
MONTHLY SUMMARY REPORT - MAY 2023



Lab Name **MEDICARE DIAGNOSTIC CENTER**

Lab No **15263**

Constituent Group **Chemistry II**

Date of Result Entered : **20/05/2023**

PT Item **Lyophilized human serum based**

Date of Report Published : **05/06/2023**

Sl No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants		Your Value	SDI	U
						CV	SD			
1	GLUCOSE II	GOD-POD II	Biosystems	310	186.96	8.32	15.56	185 mg/dL	-0.13	1.77
2	UREA II	Others (Urease Berthelot / Ned Dye) II	Any Analyser (Automation / Semi Automation)	891	25.37	18.18	4.61	26.2 mg/dL	0.18	0.31
3	CREATININE II	Jaffes Kinetic-Alkaline Picrate II	Biosystems	236	6.66	12.14	0.81	5.9 mg/dL	-0.94	0.11
4	T.BILIRUBIN II	Diazonium Salt (Colorimetric) / Jendrassik II	Biosystems	230	2.39	15.66	0.37	2.8 mg/dL	1.10	0.05
5	T-PROTEIN II	Biuret - Colorimetric II	Biosystems	278	5.78	9.57	0.55	6.8 g/dL	1.84	0.07
6	ALBUMIN II	BCG - Colorimetric II	Biosystems	268	3.38	9.52	0.32	3.8 g/dL	1.30	0.04
7	URIC ACID II	Enzymatic / Uricase Colorimetric II	Biosystems	260	6.94	11.17	0.78	7.3 mg/dL	0.46	0.10
8	CHOLESTROL II	CHOD-PAP II	Biosystems	289	121.21	9.83	11.91	169 mg/dL	4.01	1.40
9	TRIGLYCERIDE II	GPO-PAP / Enzymatic Colorimetric / End Point II	Biosystems	268	132.35	11.45	15.15	191.6 mg/dL	3.91	1.85

SDI Range	Interpretation
Within -1.00 to +1.00	Excellent.
Within ±1.01 to ±2.00	Good.
Within ±2.01 to ±2.99	Accept with caution. Warning Signal.
Beyond ±3.0	Unacceptable performance. Action Signal.

LAB ADDRESS :

MEDICARE DIAGNOSTIC CENTER
12 PVL COMPLEX, KANJIKOIL ROAD CHITTODE,
ERODE
TAMILNADU638102

Coordinator Contact Details:
Email:clinqc@cmcvellore.ac.in
Contact Number: 0416-2283102

Pamela Christudoss

Dr. Pamela Christudoss
CMC EQAS Coordinator
Christian Medical College, Vellore

Homogeneity and Stability of the sample is passed.
Data in CMC EQAS reports is confidential
CMC EQAS does not sub contract any components
***** End of Report *****