District Govt. Hospital Betul Central Processing Lab Dr. Bhimrao Ambedkar Square Betul MP 460001

Check list of Investigation/Corrective Action taken for EQAS failure

Month: 13-06-2023		
Par	ameter: Retic% (Slide)	Status
SN	Pre-analytical phase of analysis	Yes
1	Was the EQAS sample stored at the proper temperature following receipt	Yes
2	Was the test reagent stored correctly in appropriate temperature	Yes
3	Was both slides unbroken while received	100
	Analytical phase of analysis	
		Yes
1	Was the sample at room temperature	Yes
2	Was the person running the EQAS sample was trained	Yes
3	the day that the EQAS sample was run	Yes
4	Was IOC within an acceptable range on the day that the EQAS sample was run	100
•	Post-analytical phase of analysis	×
	1 of the last data of submission	Yes
1	Was the result uploaded before the last date of submission	Yes
2	Have the results been reported correctly (Match instrument raw data)	Yes
3	Was the configuration correct (instrument, method and reagent)	Yes
4	Was microscopy done again (in case of unacceptable results in Microscopic examination)	
-		

Note: Outlier result of Retic% (Slide) Z Score in parameter, result has been assesst by two slides by the pathologist, after root cause analysis has been done as per above checklist. The first result is 7.5 and the second result 7.7, and sum of two results 15.2, now the Retic% is acceptable.

डॉ. क्यू ए, बामले. पर्याताविस्ट Reviewed by कत्सालय, बे तूल

Date: 04-07-2023