

Transasia Bio-Medicals Ltd., 1st Floor, 1-10-27/1, Street No. 3, Ashok Nagar, Hyderabad - 500 020 Tel: 040-40159005/44309000 Fax: 040-44309000 Email : sz4@transasia.co.in



Brno, 10.11.2016

Laura Smart: Calibration procedure

To Whom It May Concern:

We Erba Lachema s.r.o. with registered office in Czech Republic, Karásek 1d, 621 00 Brno, represented by its Global Product Manager Pavel Suchan hereby confirm as producer of the Laura Smart instrument that this device is firmly calibrated by manufacturer and does not require any regular calibration by the user. Correct calibration of the instrument is controlled by grey strips provided initially together with the instrument.

Erba Lachema s.r.o.

Karásek 1d, 621 00 Brno DIČ: CZ26918846 -10-

RNDr. Pavel Suchan Ph.D.

Global Product Manager, Urine Analysis

Erba Lachema s.r.o.

Karásek 1d, Brno 62100, Czech Republic



