



PRIMAN INSTRUMENTS PVT. LTD.

M31112041N

Quality Product and Service for Medical Laboratory & Biotech Instruments, Consumable & Reagents "An ISO 9001:2015 Certified Company"

SAGAR GALAXY, IInd Floor, LSC, Mayur Vihar - II, New Delhi -110091
E-mail : priman@primanindia.com / email : admin@primanindia.com

Tel.: 011-41004431,
Website : www.primanindia.com

Date: 14th AUG-2023

CALIBRATION CERTIFICATE

This is to certify that Instrument model Fully Automated Biochemistry Analyzer Model **BS-430 (Sr. No. YX-14005478)** Installed at **Dr. Pankaj Diagnostics & Research Center Pvt. Ltd. Address: B- 27, Sachiwalay Colony near Kendriya Vidyalay Patna, Bihar - 800020**, has been calibrated by MINDRAY Biochemistry Calibrator MULTI SERA CALIBRATOR (Lot No.150121013 EXP:28-12-2023). During Calibration, Instrument connected through UPS.

Lamp Voltage: 12 V checked and found 11.70

Gain of all Filters are 48000- 60000.

As per the Manufacturer's Guideline the gain should be in range of
(40,000- 60,000)

Next Calibration for the Instrument is due after six month from the date of this Certificate.

For Priman Instruments Pvt. Ltd.

Authorized Signatory

Date: 14/08/2023

Place: NEW DELHI

Works :Plot No.: 75, Ground Floor, Patparganj Industrial Area, Patparganj, Delhi - 110092

ServiceUser

14-08 01:26 PM



ent Procedure

Current Step

6

check if there is water in cuvette #81-#85. If there is not, manual fill 200ul water in each cuvette; if there is, select Continue to proceed.

Gain Adjustment Result

Channel	Water Blank AD	Gain Parameters	Dark Current
340nm	49390	126	96.00
380nm	50633	61	131.00
412nm	50033	60	133.00
450nm	50311	75	119.00
505nm	49897	72	111.00
546nm	50092	79	114.00
570nm	49687	79	113.00
605nm	49527	55	138.00
660nm	50100	72	112.00
700nm	49072	75	109.00
740nm	49937	76	118.00
800nm	50124	65	118.00

Done

Continue

Stop

Activate Windows



