

S J V Scientific Company

20, Railway Lane Street, 1st Cross, Chekkan Thottam, Ganapathy, Coimbatore - 641 006.

E-mail: sjv.scientific@gmail.com. Ph no: 97866-90789

GSTIN: 33AJXPR6161K1ZG

Ref -2023/Cal/PHC/VP/013/01

To,

The Medical Officer,
Govt Primary Health Centre,
No 4, Veerapandi, Naikanur,
Coimbatore - 641 019.

CALIBRATION REPORT

MICROSCOPE CALIBRATION:

Model No : Zoom - ZBCM-08

Batch No / Serial No : ZSW/2019-20/05 / 474370

Eyepiece Lens: the lens at the top that you look through, usually 10x or 15x power.

Tube: Connects the eyepiece to the objective lenses.

Base & Arm: The bottom of the microscope, used for support. Supports the tube and connects it to the base.

Illuminator: A steady light source (110 volts) used in place of a mirror. If your microscope has a mirror, it is used to reflect light from an external light source up through the bottom of the stage.

Stage: The flat platform where you place your slides. Stage clips hold the slides in place. If your microscope has a mechanical stage, you will be able to move the slide around by turning two knobs. One moves it left and right, the other moves it up and down.


Revolving Nosepiece or Turret: This is the part that holds two or more objective lenses and can be rotated to easily change power.

Comments: This above-mentioned product is working as expected.

Calibration Date : 20.09.2023

Due Date : 19.09.2024

SJV Scientific Company



Sales & Service:

Blood Glucometer - BP Apparatus - Weighing Scale - Nebulizer - ECG Machine - Lab Kits & Chemicals - Surgical Instruments.

Sales & Service:



Scanned with OKEN Scanner