

2 & 3, Devi Illam, Priya Nagar, Telungupalayam Pirivu, Coimbatore - 641 039, TamilNadu.

T. ARUN PRASAD, M.Tech., MBA. (Ph.D)

Managing Director

Consultant Biomedical Engineer

GSTIN: 33APLPA1828A1ZM

## Calibration Protocol & Report for Medical Equipment

Equipment: CENTRIFUGE

Protect Class/Type: Risk class I

Make

LABOTECH

S.No

AV/CEN/001

Model

**BDI-152** 

Location

**ARATTAVADI - GPHC** 

S.no	Measurement Parameters			Status	Remarks
ı	Electrical Test	(b). Equip (c). Risk of (d). Chas (e). Earth	ty Standards (range: 210-240 V) potential ground & Earthing of shock/fire Hazards, if any sis & Insulation Resistance of & Chassis Leakage Current ry condition & charging status	230volts 0 to 5volts No Zero ohms 10mA NA	Stabilizer/UPS With in ranges well
H	Physical Condition Test	(a). Cabir (b). Front asser	(a). Cabinet & Visible damages, if any (b). Front panel/top lid, buckets, motor assembly & Power cord (c). Risk of mechanical damages, if any		
111	Performance Test	(a). On/o (b).Speed functi (c). Runni (d).Abnor (e). Displa (f). Keypa	<ul> <li>(a). On/off switch</li> <li>(b). Speed regulator switch &amp; functional tests</li> <li>(c). Running test</li> <li>(d). Abnormalities</li> <li>(e). Display checkup</li> <li>(f). Keypad/indicators &amp; connectors</li> </ul>		ok perfect perfect
IV	Calibration Tes	(b). Speed wobbl (c). Runni (d). Displa	<ul> <li>(a). Vibration test</li> <li>(b). Speed regulation gradual hike/ wobbling tests</li> <li>(c). Running test with load</li> <li>(d). Display checkup</li> <li>(f). timer test</li> </ul>		perfect 3600 RPM
Formulated By: Einstein Research Laboratory			issued By:	User Department:	
Endorsed By: Research team head  Reviewed On: 12.09.23 Next Review: 12.03.24		3	T.Arunprosad SILICON TECHNOLOGIES A Biomedical Instrumentation Compan COIMBATORE - 641 039		

Cell: 98652 85754, 81246 75020.

Email: st.biomedical@gmail.com

**5** 98652 85754, 81246 75020