

MAKE

G MULTI SERVICES

78, Yahoob Street, Gobichettipalayam, Erode Dist, Tamil Nadu. Pin: 638 452. gmultiservicesgobi@gmail.com. Cell: 96597 20654. **Bio Medical Instruments Services**

Er. A. Ibathul Raziq DECE., BE (ICE) **Bio Medical Engineer**

Calibration Protocol & Report for Medical Equipment

Equipment: CELL COUNTER

: SYSMEX

: XP-100 MODEL

PROTECT CLASS / TYPE : Risk Class I

: B5885

S. No. Location

: UGPHC, KUGALUR.

S.No	Me	asurement Parameters	Status	Remarks
l	Electrical Test	(a)Safety Standards (range: 210-240 V) (b)Equipotential ground & Earthing (c)Risk of shock/fire Hazards,if any (d)Chassis & Insulation Resistance (e) Earth & Chassis Leakage Current (f) Battery conditions & charging status	230 volts 0 to 5 volts No Zero ohms 10mA NA	Stabilizer/UPS With in Ranges Well
2	Physical Condition Test	(a)Cabinet & Visible damages, if any (b) Front Panet/Top lid, buckets,motor aseembly & power cord (c)Risk of mechanical damages,if any	No OK No	Checked OK
3	Performance Test	(a) On /Off switch (b) switch & functional tests (c) Spindle Motor (d) Running test (e) Abnormalities (f) Display checkup (g) Keypad/indicators & connector	Normal OK OK OK No Well Well	Ok Perfect Perfect
4	Calibration Test	(a) WB Test (b) Printer (c) Display checkup (d) Timer test (e) Result	Done Ok OK OK Perfect	Perfect
Formulated By: Edition Research Laboratory Endorsed By: Research team head Reviewed On: 18.08.23 Next Review: 17.08.24		Issued By:	User Depa	rtment: