

• Remove the cap
• Remove the cap and use before
• Cover the spill with absorbent cotton or tissue paper
• Discard this tissue paper in Yellow Bag
• Disinfect the surface with 1% Sodium Hypochlorite solution, 8% Phenol, 10% Bleach for 10-15 minutes or use Phenolic Disinfectant
• Now wash the surface with soap and water
• Remove PPE and discard in waste yellow bag

SUNSHINE BIOMEDICAL SOLUTIONS
Mail: ldcalibration@sunshinebiomedical.com
Mob: +91 73587 77134
Calibration

Equipment Name: Centrifuge
Serial No: 20111925
Date: 06-10-23 Due: 05-10-24
By: 14.11.23

ISO 9001 : 2015