

TIN : 07650135381

NEW GOLDEN ENTERPRISES

STOCKISTS



BDH, Glaxo, S.D. Fine, Merck, Finar, Borosil Glassware, Whatman Filter Paper
Tarson/Luxbro, Diagnostic Kits Imported Chemicals & General order Suppliers
Office : Post Box No. 2177, Delhi-110007

8 UA, Ground Floor, Jawahar Nagar, Near Kamla Nagar, Delhi-110007 (India)
E-mail : newgolden_vinay@rediffmail.com




Certificate of Calibration (Model – Mindray BC3600)

Date – 05.03.2023

This is to certify that the below mentioned instrument in Patliputra Nursing Home, Samaypur - Libaspur Rd, Samaipur, Yadav Nagar, Samaypur, Delhi, 110042 has been calibrated by us 05th March'2023. The next calibration is due on before 04th March'2024. The values are falling within the specified range.

- Instrument Model No. BC3600
- Model Name – MINDRAY BC3600
- Lot NO. PLUS0423
- Exp. : 05-05-2023
- Sr. No. TB-89002330

For New Golden Enterprises

 (Service Engineer)

SC-CAL PLUS

HEMATOLOGY CALIBRATOR

CALIBRATOR
LOT PLUS0423

2023-05-05

Instrument	Parameter	Calibration Values	Acceptable Range
BC-3600, BC-3300	WBC x 10 ⁹ /L	9.4	9.2 - 9.6
	RBCx 10 ¹² /L	4.53	4.45 - 4.61
	HGB g/L	130	128 - 132
	MCV fL	84.6	82.6 - 86.6
	PLT x 10 ⁹ /L	265	253 - 277
BC-3200, BC-3000 CT	WBC x 10 ⁹ /L	9.7	9.5 - 9.9
	RBCx 10 ¹² /L	4.46	4.38 - 4.54
	HGB g/L	130	128 - 132
	MCV fL	85.9	83.9 - 87.9
	PLT x 10 ⁹ /L	256	244 - 268
BC-3000 Plus, BC-2900 BC-1800	WBC x 10 ⁹ /L	9.4	9.2 - 9.6
	RBCx 10 ¹² /L	4.57	4.49 - 4.65
	HGB g/L	133	131 - 135
	MCV fL	84.6	82.6 - 86.6
	PLT x 10 ⁹ /L	249	237 - 261
BC-2800, BC-2600 BC-2800Vet, BC-2600Vet	WBC x 10 ⁹ /L	9.1	8.9 - 9.3
	RBCx 10 ¹² /L	4.43	4.35 - 4.51
	HGB g/L	131	129 - 133
	MCV fL	86.2	84.2 - 88.2
	PLT x 10 ⁹ /L	252	240 - 264
BC-2300, BC-2100	WBC x 10 ⁹ /L	9.0	8.8 - 9.2
	RBCx 10 ¹² /L	4.35	4.25 - 4.45
	HGB g/L	129	127 - 131
	MCV fL	85.5	83.5 - 87.5
	PLT x 10 ⁹ /L	245	230 - 260
BC-20s, BC-21s BC-30, BC-30s, BC-30e BC-31, BC-31s	WBC x 10 ⁹ /L	9.2	9.0 - 9.4
	RBCx 10 ¹² /L	4.33	4.23 - 4.43
	HGB g/L	131	129 - 133
	MCV fL	87.3	85.3 - 89.3
	PLT x 10 ⁹ /L	249	234 - 264
BC-10, BC-10e, BC-11 BC-20, BC-21	WBC x 10 ⁹ /L	9.0	8.8 - 9.2
	RBCx 10 ¹² /L	4.34	4.24 - 4.44
	HGB g/L	130	128 - 132
	MCV fL	87.6	85.6 - 89.6
	PLT x 10 ⁹ /L	245	230 - 260

Before using, refer to the instruction sheet for mixing directions.


Mindray

Shenzhen Mindray Bio-Medical Electronics Co., Ltd.

Mindray Building, Keji 12th Road South, Hi-tech Industrial Park, Nanshan, ShenZhen 518057, P.R.China

Tel: +86 755 81888998

Fax: +86 755 26582680

EC REP

Shanghai International Holding Corp. GmbH (Europe)

Eiffestraße 80 20537 Hamburg, Germany

Tel: 0049-40-2513175

Fax: 0049-40-255726

Installation Certificate For MINDRAY(BC-3600) Sr. No.- TB-89002330

This is to certify that the instrument MINDRAY(BC-3600) Sr. No.- TB-89002330 is successfully installed Patliputra Nursing Home Samaypur - Libaspur Rd, Samaipur, Yadav Nagar, Samaypur, Delhi, 110042 and comminssioned at :and the Installation Protocol/ checklist has been Successfully completed for the above Instrument.

Date of Installation : 07.05.2018

Patliputra Nursing Home

Samaypur - Libaspur Rd, Samaipur,
Yadav Nagar, Samaypur, Delhi, 110042

Name : Shivam Gupta

Designation : Service Engineer

Installation Qualification for MINDRAY(BC-3600) Sr. No.- TB-89002330

**Customer Name : Patliputra Nursing Home Samaypur - Libaspur Rd, Samaipur,
Yadav Nagar, Samaypur, Delhi, 110042**

Instrument Name : MINDRAY(BC-3600)

Serial No : TB-89002330

Initial Inspection of the unit carried out and the details are as follows:

Shivam



System Condition Report:

Found the system to have been delivered in satisfactory condition and no external physical damaged observed on the same, Package was kept in good Condition as per the directional indicators like not tilt, indicating the system has not been subjected to mechanical shocks or stored in any manner, so as to cause any damage to the same.

Found all the required accessories as required.

Intallation Procedure And Checklist Attached for the records.

External Requirement for Installation:

1. Input voltage of 220V-240V/50Hz or 60 Hz.
2. Recommended operating Temperature is 15-30 degree Celsius, with in Relative Humidity 30-85% and Atmospheric pressure 70-106kPa.

Installation Qualification for MINDRAY(BC-3600) Sr. No.- TB-89002330

Carried out all the installation procedures as per the installation procedure and checklists.

Carried out all the necessary checks and alignments.

Carried out all the necessary system checks and tests.

Handover the Instrument for operatos Training And Qualifications



For New Golden Enterprises

A handwritten signature in blue ink that reads "Shivam".

Name : Shivam Gupta

Designation : Service Engineer

Performance Qualification for MINDRAY(BC-3600) Serial No: TB-89002330

Calibration Parameters

Checked and found all the precision check also.

System Certification:

Study data has determined that the system described in this document either meets the necessary criteria outlined in this Performance Qualification Protocol, or exceptional conditions have been identified and documentation included.

The system is ready for Specific Usage.

Protocol Performed By:

Name : Shivam Gupta

Designation : Application Manager

PERFORMANCE QUALIFICATION DONE BY

Name : Shivam Gupta

Designation : Application Manager



Al. vany

Operational Qualification:

System Certification:

Study data has determined that the system described in this document either meets all criteria outlined in this Operational Protocol, or exceptional conditions have been identified and documentation included.

Exceptional Conditions, If any Have been Addressed.

The System is ready for specific usage.

Protocol Performed by : **New Golden Enterprises**



Name: **Shivam Gupta**

A handwritten signature in blue ink that reads "Shivam".

Designation: Service Engineer

Customer Authorization Patliputra Nursing Home Samaypur
Libaspur Rd, Samaipur, Yadav Nagar,
Samaypur, Delhi, 110042

Name: //

Designation: **Manager**

Company Representative Name And Sign

//
Customer Name And Sign

Date:

Date:

Operations Qualifications for MINDRAY(BC-3600) Sr. No.: TB-89002330

1. Verified all the Mechanical Movements : Done
2. Verifeied Hydraulics : Done
3. Verified Electrical Systems : Done
4. Verified the all reagents Systems : Done

NEW GOLDEN ENTERPRISES

Name : Shivam Gupta

Designation : Service Engineer

Date : 07.03.2023 Instrument Setup

1. Assembled the instrument accessories.
2. Removed the shipping Clamps.
3. Connected the LYSE
4. Connected the RINSE
5. Connected the DILUENT REAGENT.
6. Connected the Waste Tubing.
7. Connected the Power cord and connection cord.
8. Coonected the External Printer.
9. Install the software.
10. Initialise the machine and follow the installation procedure.



Operational Instection

1. Checked and found Mechanical Movements OK.
2. Checked and found Hydraulics OK.
3. Checked and found Electricals OK.
4. Checked with Controls And Samples, Results Are found OK.



New Golden Enterprises

Name : Shivam Gupta

Designation : Service Engineer