

78, Yahoob Street, Gobichettipalayam, Erode Diet, Tamil Nadu. Pin:638 452. gmultiservicesgublingmail.com. Cell: 96597 20654. Bio Medical Instruments Services

Er. A. Ibethul Raziq DECE., BE (ICE) Bio Medical Engineer

Calibration Protocol & Report for Medical Equipment

Equipment : CENTRIFUGE MAKE : NAUTH

S. No.

PROTECT CLASS / TYPE : Risk Class !

MODEL : CENT 100 Location

: AR/CENT/001 : GPHC, ARACHALUR

S.No	Measurement Parameters		Status	Remarks
1	Electrical Test	(a)Safety Standards (range: 210-240 V) (b)Equipotential ground & Earthing (c)Risk of shock/fire Hazards, if any (d)Chassis & Insulation Resistance (e) Earth & Chassis Leakage Current (f) Battery conditions & charging status	230 volts 0 to 5 volts No Zero ohms 10mA NA	Stabilizer/UPS With in Ranges Well
2	Physical Condition Test	(a)Cabinet & Visible damages, if any (b) Front Panet/Top lid, buckets,motor aseembly & power cord (c)Risk of mechanical damages, if any	No OK No	Checked DK
3	Performance Test	(a) On /Off switch (b) Speed regulator switch & functional tests (c) Running test (d) Abnormalities (e) Display checkup (f) Keypad/indicators & connector	Normal OK Accurate OK No Well Well	Ok Perfect Perfect
4	Calibration Test	(a) Vibration test (b) Speed regulation gradual hike/wobbling test (c) Running test with load (d) Display checkup (e) Timer test	Done Ok	Perfect
Formulated By Edition Research Laboratory		Issued By:	User Department:	
Resear head Review	sed By : rch team wed On: 02.10.23 leview : 01.10.24	G. A. JEJ		

The Chille of the Constraint of the Chille of the Constraint of the Constraint of the Chille of the