

AUGUST
23

Glucose	194	2.05	
Creatinine	1.9	1.18	
Cholesterol	124	2.20	

CAPM

- * Creatinine - SDI Range within limit
- * Glucose and Cholesterol - SDI Range accept with caution - warning signal
- * Called BME - maintain SAA cleaning tube and washing also check Standard
- * Recalibrated Glucose & Cholesterol test
 Rerun PQAS and EQAS

IQAS - Glucose value	- 106
IQAS - Cholesterol value	- 110
EQAS - Glucose value	- 178
EQAS - Cholesterol value	- 108