



Metropolis Healthcare Ltd., EQAS

PRIMUS PATHOLOGY & DIAGNOSTICS (1015)

Outlier And Analyte Summary Report

Outlier Details For Cycle No 230310 and Sample No 01

Report Date : 19/12/2023

Analyte	Instrument	Result Value	Standard Unit
X Blood (urine)	Manual	4+	-
Legend @ : Acceptable X : Unacceptable # : Not Evaluated ⌚ : Delayed Result Entry			

Problem Classification: _____

Corrective Action: _____

Reviewed by: _____


Dated: _____

Instrument : Manual

Analyte	Standard Unit	Result Value	Accepted Value
@ pH	-	6.5	6,6.5,7
@ Specific Gravity	-	1.015	1.015,1.020,1.025
@ Glucose (urine)	-	3+	+2,+3,+4
@ Protein (Urine)	-	1+	+1,+2,+3
@ Ketone (urine)	-	Negative	Neg
X Blood (urine)	-	4+	+1,+2,+3
@ Bilirubin (urine)	-	Negative	Neg
@ Urobilinogen (urine)	-	Normal	Nor

Legend @ : Acceptable
 @ : Acceptable
 X : Unacceptable
 # : Not Evaluated
 Ⓞ : Delayed Result Entry

Total Parameters	8
Not Evaluated Parameters	0
Evaluated Parameters	8
Outlier Parameters (X)	1
EQAS Score Clinical Path - Urine Routine Chemistry	87.50 %



Dr Puneet Kumar Nigam
 PT coordinator & Technical Manager, MHL EQAS
 Unit No. 409-416.
 Commercial Building - 1A
 Kohinoor Mall, Kiroil Road, Kurla (W),
 Mumbai - 400070

CORRECTIVE & PREVENTIVE ACTIONS

ROOT CAUSE ANALYSIS, CORRECTIVE AND PREVENTIVE ACTION (CAPA)

Non Compliance Reference

No. _____ dated _____

Dept: _____

Date: _____

Incident	EQAS - Clinical Pathology - Urine Blood was unacceptable.
Root cause analysis/Investigation	Urine blood finding by strip showed 3+ but slide interpretation done by the staff was 4+.
Tool Used/ Risk Assessment	-
Immediate Correction	-
Corrective Action	Retraining of all senior staff on slide interpretation for urine was done. Spot sample analysis was done by Dr. Neeraj for all senior staff & reporting was compared.
Preventive Action	Slide interpretation will be crosschecked in such cases.
Verification of CAPA Done	Done By: _____ Date & Sign: _____

Date: 23/12/2023

Subject : Slide analysis Comparison of Senior Staff

Patient Name: Aarti Pandey

Test : Microscopy results of Urine Sample

Urine Sample Parameters	Dr. Neeraj	Ms. Sushma	Ms. Tabassum
Pus Cell	4-6	3-4	
Epithelial Cell	3-4	2-3	
RBC	Absent	Absent	
Other	Yeast cells seen	Yeast organism	

Comment: Slide analysis comparison of staff has noted to be acceptable.