



CERTIFICATE OF CALIBRATION

ADL/BO/MS/C/001/24-25

INSTRUMENT	BIOLIS 30I PREMIUM Fully Automated Clinical					
Date of Installation: 18/07/2022	SITE : South Zone Hospital	& Researc	h Center	CITY: Kunnamkulam		
State : Kerala	Country: INDIA S/N:6010520				520422	
Section	Nature of work	Clean / Check		Lubricate	Replace	
Reaction Tray	Cleaning of Reaction Tray & Checking of Movement	V	V			
Sampling Pump	Cleaning	V V		٧		
Reagent Pump	Cleaning	v v		V		
Probe Washing Pump	Cleaning	v v				
Reaction Cuvette Washing Pump	Cleaning	v v				
Sample Probe	Cleaning	V V				
Reagent Probe	Cleaning	v v				
	Movement Check					
Barcode Scanner	Scanning Check	٧				
Sampling Pump	Movement Check	V				
Reagent Pump	Movement Check	V				
Probe Washing Pump	Movement Check	٧				
Reaction cuvette Washing Pump	Movement Check	٧				
Sampling Probe	Movement Check	٧				
Reagent Probe	Movement Check	٧				
Mixing Unit	Up & Down Movement	٧				
Mixing Pump	Movement Check	√				
Water Supply & Drainage	Movement Check	٧				
Control PCB	Checking of Connections	٧				
Power Supply (control board)	Checking O/P Voltage(5.26 V)	٧				
Body	Checking of Leakage	٧				
Total Inspection	Cell Blank Test	√				
Temperature Check	Reaction Tray 37 Celsius	V				
	Washing water 38 Celsius	V				

AGAPPE DIAGNOSTICS LTD. 150 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN: U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT Agappe Hills, Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India. Tel: + 91 484 286 7000 | Email: agappe@agappe.in

MUMBAI [REGISTERED OFFICE] MUMBAI (REGISTERE DEFICE) 401 & 402, 4" Floor, Jaisingh Business Centre, 119, Sahar Road, Parsiwada, Andheri (East), Mumbai - 400 099, India. Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in **EQUIPMENT PLANT** X/588-CB, Block No. 32, KINFRA Small Industrial Park, Nellad, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

DELHI OFFICE DSM 540, 5th Floor, DLF Tower, Shivaji Marg, New Delhi - 110 015, India. Tel: +91 I1 4558 8416 | Email: delhioffice@agappe.in

KOLKATA OFFICE 406, Merlin Matrix, Plot No-10, Block-DN, Sector V, Salt Lake City, Kolkata - 700 091. Tel: +91 33 4003 0451 | Email: kolkattaoffice@agappe.in

BANGALUKE OFFICE 5-6, II* Floor, Red Cross Bhavan, No:26 Race Course Road, Bangalore - 560 001. Tel: +91 80 2228 8288 Email: bangaloreoffice@agappe.in





	Wave Length	OD Data	Check
Auto Gain Check	340nm	49239	√ OK
	380nm	49169	√ OK
	405nm	48953	√ OK
	450nm	50997	√ OK
	505nm	48997	√ OK
	546nm	49261	√ OK
	570nm	51054	√ OK
	600nm	51029	√ OK
	660nm	49252	√ OK
	700nm	50725	√ OK
	750nm	50789	√ OK
	800nm	48981	√ OK

The equipment is calibrated for the different parameters by using Multi Calibrator

This is to certify that the above-mentioned product has been calibrated according to the standard operating procedure given by the Tokyo Boeki Medical Systems, Japan.

Date of Calibration: 20th Feb 2024 Next Due on: 19th July 2024

For AGAPPE DIAGNOSTICS LTD.,



Authorized signatory

Name: Mahmood Basheer P.P.

Designation: Product Specialist, Mispacare

AGAPPE DIAGNOSTICS LTD. ISO 9001:2015 | EN ISO 13485:2016 CERTIFIED COMPANY | CIN: U24239MH1998PLC115413

CORPORATE OFFICE / REAGENT PLANT Agappe Hills, Pattimattom (PO), Dist. Ernakulam, Kerala - 683 562, India. Tel: + 91 484 286 7000 | Email: agappe@agappe.in

MUMBAI [REGISTERED OFFICE]
401 & 402, 4" Floor, Jaisingh Business Centre, 119,
Sahar Road, Parsiwada, Andheri (East), Mumbai - 400 099, India.
Tel: +91 22 4300 8000 | Email: mumbaioffice@agappe.in

EQUIPMENT PLANT X/588-CB, Block No. 32, KINFRA Small Industrial Park, Nellad, Cochin, Kerala, India - 686 721. Tel: +91 484 276 7477.

DELHI OFFICE DSM 540, 5" Floor, DLF Tower, Shivaji Marg, New Delhi - 110 015, India. Tel: +91 11 4558 8416 | Email: delhioffice@agappe.in KOLKATA OFFICE 406, Merlin Matrix, Plot No-10, Block-DN, Sector V, Salt Lake City, Kolkata - 700 091. Tel: +91 33 4003 0451 | Email: kolkattaoffice@agappe.in

BANGALORE OFFICE
5-6, II* Floor, Red Cross Bhavan, No:26 Race Course Road,
Bangalore - 560 001. Tel: +91 80 2228 8288
Email: bangaloreoffice@agappe.in