

FUJIFILM

Value from Innovation

Date:03.07.2023

To,
Dr. Suresh Teckchandani
Vinamra Swaraj Hospital, Plot No. 108/109, Opposite Shiv Sena Shaka ,Sector 11, Vashi, Navi Mumbai - 400 703.Navi
Mumbai - 400 703.

Dear Sir/Madam

We thank you for your trust towards Fuji Medical Products and Service Delivery. Hope you have had a good experience with the product and service. In order to continue our best service delivery, we propose the service contract for following equipment's installed at your centre/hospital are as follows.

Equipment	FDC NX 700i
Serial No	97021446
Installation Date	13 Jul 2020
Existing Contract period end date	12 Jul 2023
Customer Code	1012792
Type of Contract	Annual Maintenance Contract
New Contract Period Start-End Date	13 Jul 2023 to 12 Jul 2024
Basic Price Rs.	30000/-
GST Rs. @ 18%	5400/-
Net Value Rs.	35400/-

Terms and Conditions: AMC

1. 100 % advance payment along in favor of FUJIFILM INDIA PRIVATE LIMITED -- RTGS detail- HDFC Bank, Vatika Atrium, A-block, Golf Course Road, Sector-53, Gurugram – 122002, IFSC Code of the Branch - **HDFC0000572** (HDFC Bank), MICR CODE-**110240093**. Account Number **FUJI401012792**
2. 100% payment on basis of **Proforma invoice** and PM report submission.
3. Any failure due to Third party equipment like UPS, batteries or earthing will not be considered as down time of machine.
4. Contract includes 4 preventive maintenance calls. Any number breakdown calls are covered
5. No Labour charges, whereas spares will be on chargeable basis
6. All the Defective replaced spare parts and Good Spare parts to be return back to Fuji warehouse
7. For services under this contract customer can contact Toll Free Number given below or FFIN Engg. Mr. Prabhakar Sharma - 9699258155

Thank You



Ritvik Sen
Regional Service Manager West
Fujifilm India Pvt. Ltd.

Dr. S.D. Teckchandani

Consulting Pathologist

VINAMRA SWARAJ HOSPITAL
Sec:11,Vashi, Navi Mumbai - 400 703

Toll free assistance line is now open to you from Mon - Sat (09:00 to 05:30 pm), for better service please log all service complaints @ 1800 103 3854 or email us at medical.services@fujifilm.com

FUJIFILM INDIA PRIVATE LIMITED

Corporate Office : Unitech Cyber Park, 8th floor, Unit no.801-807, Tower C, Sector-39, Gurugram- 122101, Haryana Telephone : +91-124-4325500, Fax : +91-124-4325555, Website : www.fujifilm.com , Email : Contact@fujifilm.com

Registered Office : Business Centre Office No. 521, Level V, Caddie Commercial Tower, Hospitality District Aerocity IGI Airport, New Delhi-110037, Telephone : +91-11-66544932.

Mumbai Office: Unit No.504 & 505, 5thFloor,349 Business Point, Western Express Highway, Andheri (East), Mumbai-400069
Telephone: +91-22-42364000, Fax +91-22-42364001.

Subject : No Equipment Calibration Requirement for FDC 4000i, NX500i, FDC7000i and NX700i

Design Change Reason	Action	Time frame
<input type="checkbox"/> Specification/Function change <input type="checkbox"/> Pro-active improvement <input type="checkbox"/> Serviceability <input type="checkbox"/> Part Replacement <input type="checkbox"/> Corrective Action <input type="checkbox"/> Preventive Action <input type="checkbox"/> Regulatory Compliance <input checked="" type="checkbox"/> Other (information)	<input type="checkbox"/> Mandatory <input type="checkbox"/> Recommended <input type="checkbox"/> Not required () <input checked="" type="checkbox"/> Other (information only)	<input type="checkbox"/> Immediately <input type="checkbox"/> Next maintenance <input type="checkbox"/> Upon trouble situation <input type="checkbox"/> At your convenience <input checked="" type="checkbox"/> Not Applicable <input type="checkbox"/> Other ()

Date: 04/December/2014 Dept.: Medical Systems Business Div.

Originator: M. Jumawid Edited by: M. Jumawid Approved by: Y. Kashi

Affected Products: FUJI DRI-CHEM analyzers FDC 4000i, NX500i, FDC7000i and NX700i

To whom it may concern:

This document certifies that FUJI DRI-CHEM analyzers FDC4000i, NX500i, FDC7000i and NX700i do not require equipment calibration.

The quality of these analyzers is guaranteed through the following principle.

FUJI DRI-CHEM analyzers are equipped with reference white and black plates. The functions of these plates are 1) to adjust the gain for all the interference filters before every measurement, and 2) to assure that the analyzer is working within the manufacturer's quality standard.

Measurements are conducted after gain adjustment, where the optical densities of reflectance for these plates measured in a constant interval, are checked. If the values exceed the threshold, the machine will indicate an error and the analyzer will not give a measurement result for this abnormal case.

The reference white and black plates are made from ceramic material which do not degrade over time, thus the quality of the analyzers can be confirmed over a long period of time.

- Special tools / Instruments : Yes No
- Supply of the articles : Yes No
 If Yes, when :
 Name / Code :
- Old parts after replaced : Yes No
- Report requirement after finished : Yes No

Contact person: Mariejoy Jumawid Tel: +81-3-6418-2213 E-mail: mariejoy.jumawid@fujifilm.com