

CALIBRATION CERTIFICATE

To,

Welldynamic Diagnostic - Rohini Delhi

Instruments : AutoQuant-100 Amara Machine Serial Number : **AQ101073**

Issued On: **10-Apr-2024** Valid Upto: **09-Apr-2025**

Appearance	No Imbalance Probe Movement
Leakage in Water tubing or CRU	No water leakage
Ambient temp (15-32) Degree	23 Degree
Lamp Check	Ok
Threshold Value	34,500 (≥25,000)
Sample Dispersing Volume	In range (0 to 30 ul)
Reagent cooling Temperature	7.5 Degree (0 to 15 Degree)
Optical Path length	0.65
Waster & Water Sensor	Working Properly
Valve & Motor	Working Properly
Reagent Dispersing Volume	In range (0 to 150 ul)
Power Supply and Earthing	230 Volts and Grounding less then 5 V

Conclusion: All parameter are within range. The Machine is working Properly. (As per Meril terms.)

Engineer Name: Mr. Devindra Kumar

Engineer Sign:

Stamp:



T +91 260 3052100 F +91 260 3052125

www.merillife.com

Meril Diagnostics Pvt. Ltd. Survey No.135/2/B, Bilakhia House Muktanand Marg, Chala Vapi - 396191. Gujarat, India.

Cardiovascular Orthopedic Surgical Diagnostics