
	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

**Section 1 - Initiation of ANF (to be filled by the person who raising the ANF)**


PT/EQAS Agency <input type="checkbox"/> CAP <input type="checkbox"/> AIIMS <input type="checkbox"/> BIORAD <input type="checkbox"/> CMC <input type="checkbox"/> RML <input type="checkbox"/> Other				
Survey name & distribution ID: Metropolis helathcare LIMITED EQAS Date on agency report: 10-July-2024				
ANF No: ANG/JUL/24/01	Issued by: Nicky Anand	Issue Date: 18-JUL-2024	Due date: 28-JULY-2024	
Department: Clinical Pathology				
ANALYTE or EXAMINATION: Specific Gravity				
Sample ID	Result submitted	PT targets	PT acceptable range	Problem/Performance
01	1.020	NA	NA	Unacceptable
Comment /Observations: (e.g. trend, previously missed within last 12 months, lab in regulatory jeopardy for this analyte?)				

**Section 2 - Investigation of Non Conformance- Checklist**

Sl	ANF-CHECKLIST	Yes	No	N/A
1	Specimen temperature check, as per kit instructions	Yes		
2	Specimen storage condition check, as per kit instructions	Yes		

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

<b>3</b>	Specimen physical condition check	<b>Yes</b>		
<b>4</b>	Sample integrity up to arrival in lab: any shipping, delay or other sample problems?		<b>No</b>	
<b>5</b>	Sample integrity in-house: any problems with sample handling in the lab?		<b>No</b>	
<b>6</b>	Were there any instrument problem?		<b>No</b>	
<b>7</b>	Were there any method problem?		<b>No</b>	
<b>8</b>	Were there any faulty reagent/QC and Calibrator?		<b>No</b>	
<b>9</b>	Were there QC trends / problems at time of assay?		<b>No</b>	
<b>10</b>	Was Peer group data checked, if required			<b>N/A</b>
<b>11</b>	Were there any Calibration (Intercept/slope) problems at time of assay?		<b>No</b>	
<b>12</b>	Was water quality checked?			<b>N/A</b>
<b>13</b>	Did any technical errors occur due to pipetting error		<b>No</b>	
<b>14</b>	Did any technical errors occur due to sample mix-up		<b>No</b>	
<b>15</b>	Did any technical errors occur due to incorrect process, other than reconstitution, dilution or calculation errors,		<b>No</b>	
<b>16</b>	Did any technical errors occur due to misinterpretations		<b>No</b>	
<b>17</b>	Was the instrument checked for daily. Weekly. Monthly, semiannually and Annual maintenance	<b>Yes</b>		
<b>18</b>	Was there a lab clerical error? (e.g. error in unit conversion needed for survey only, transcription error onto result form, wrong method details submitted, Factor removal)		No	
<b>19</b>	Was there a clerical error by the PT / EQA agency?		No	
<b>20</b>	Were there Patient data trends / problems at time of assay?		No	
<b>21</b>	Were there any gaps / issues in training or competency assessment?		No	
<b>22</b>	Was the sample condition appropriate at time of retesting (mention temperature...)		No	
<b>23</b>	Has sample been re-tested / re-examined?		No	
<b>24</b>	Has the lab perform Interlaboratory Comparison		No	

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

Questions 4-22 give details for any Yes answers

Instrument daily, weekly and monthly maintenance was done.


Question 23/24: if answer is yes, give results of repeat testing

Sample ID	Original result	Repeat value/ Lab Result	PT Targets/ Referral Lab Result	PT/ILC Acceptable Range	Status (Acceptable/Not acceptable)
<u>D5397330</u>	-	<u><math>\geq 1.030</math></u>	<u><math>\geq 1.030</math></u>	$\pm 10\%$	<u>Acceptable</u>
<u>D5397894</u>	-	<u>1.010</u>	<u>1.010</u>	$\pm 10\%$	<u>Acceptable</u>

**Section 3 Root cause (Refer to Non-conformance error Reason)**


Random error, IQC of the day satisfactory.

Why do you think the non-conformance / error occurred? Use this area to explain your findings.

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

<p>Assessment of the impact of root cause on patient results (Patient results may be affected before, during or after the PT event). Describe how the conclusion of impact was made and any corrective actions made to the patient result(s); "If No" - Explain why there was no patient impact?</p>
<p><b>IQC Values are within acceptable range.</b></p>
<p>Describe any previous proficiency testing issues with this test in the last sample:</p>

**Section 4: Department - Conclusion & Corrective / Preventive Action**

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>




<b>What corrective action have you carried out?</b>
Split Testing has been done , results are satisfactory.
<b>Preventive action put into place?</b>
NA
<b>Due date for closure of proposed corrective and preventive action:</b>
28-July-2024

**Person investigated**

**Department Manager**


**Lab Head**


**Signature with Date:**

Note: After signatures please hand over the form along with supporting documents to QA.

**Date when ANF received to QA along with supporting documents:** \_\_\_\_\_ **by:** \_\_\_\_\_

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>


	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

**Section 1 - Initiation of ANF (to be filled by the person who raising the ANF)**

<u>PT/EQAS Agency</u> <input type="checkbox"/> CAP <input type="checkbox"/> AIIMS <input type="checkbox"/> BIORAD <input type="checkbox"/> CMC <input type="checkbox"/> RML <input type="checkbox"/> Other				
<u>Survey name &amp; distribution ID:</u> Metropolis helathcare LIMITED EQAS <u>Date on agency report:</u> 10-July-2024				
<u>ANF No:</u> ANG/JUL/24/02	<u>Issued by:</u> Nicky Anand	<u>Issue Date:</u> 18-JUL-2024	<u>Due date:</u> 28-JULY-2024	
<u>Department:</u> Clinical Pathology				
<u>ANALYTE or EXAMINATION:</u> PH				
Sample ID	Result submitted	PT targets	PT acceptable range	Problem/Performance
01	5.5	NA	NA	Unacceptable
<u>Comment /Observations:</u> (e.g. trend, previously missed within last 12 months, lab in regulatory jeopardy for this analyte?)				


**Section 2 - Investigation of Non Conformance- Checklist**

Sl	ANF-CHECKLIST	Yes	No	N/A
1	Specimen temperature check, as per kit instructions	Yes		

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

<b>2</b>	Specimen storage condition check, as per kit instructions	<b>Yes</b>		
<b>3</b>	Specimen physical condition check	<b>Yes</b>		
<b>4</b>	Sample integrity up to arrival in lab: any shipping, delay or other sample problems?		<b>No</b>	
<b>5</b>	Sample integrity in-house: any problems with sample handling in the lab?		<b>No</b>	
<b>6</b>	Were there any instrument problem?		<b>No</b>	
<b>7</b>	Were there any method problem?		<b>No</b>	
<b>8</b>	Were there any faulty reagent/QC and Calibrator?		<b>No</b>	
<b>9</b>	Were there QC trends / problems at time of assay?		<b>No</b>	
<b>10</b>	Was Peer group data checked, if required			<b>N/A</b>
<b>11</b>	Were there any Calibration (Intercept/slope) problems at time of assay?		<b>No</b>	
<b>12</b>	Was water quality checked?			<b>N/A</b>
<b>13</b>	Did any technical errors occur due to pipetting error		<b>No</b>	
<b>14</b>	Did any technical errors occur due to sample mix-up		<b>No</b>	
<b>15</b>	Did any technical errors occur due to incorrect process, other than reconstitution, dilution or calculation errors,		<b>No</b>	
<b>16</b>	Did any technical errors occur due to misinterpretations		<b>No</b>	
<b>17</b>	Was the instrument checked for daily. Weekly. Monthly, semiannually and Annual maintenance	<b>Yes</b>		
<b>18</b>	Was there a lab clerical error? (e.g. error in unit conversion needed for survey only, transcription error onto result form, wrong method details submitted, Factor removal)		<b>No</b>	
<b>19</b>	Was there a clerical error by the PT / EQA agency?		<b>No</b>	
<b>20</b>	Were there Patient data trends / problems at time of assay?		<b>No</b>	
<b>21</b>	Were there any gaps / issues in training or competency assessment?		<b>No</b>	
<b>22</b>	Was the sample condition appropriate at time of retesting (mention temperature...)		<b>No</b>	
<b>23</b>	Has sample been re-tested / re-examined?		<b>No</b>	
<b>24</b>	Has the lab perform Interlaboratory Comparison		<b>No</b>	



	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

Questions 4-22 give details for any Yes answers

Instrument daily, weekly and monthly maintenance was done.


Question 23/24: if answer is yes, give results of repeat testing

Sample ID	Original result	Repeat value/ Lab Result	PT Targets/ Referral Lab Result	PT/ILC Acceptable Range	Status (Acceptable/Not acceptable)
<u>D5397330</u>	-	<u>5.5</u>	<u>5.5</u>	±10%	<u>Acceptable</u>
<u>D5397894</u>	-	<u>6.0</u>	<u>6.0</u>	±10%	<u>Acceptable</u>

**Section 3 Root cause (Refer to Non-conformance error Reason)**


Random error, IQC of the day satisfactory.

Why do you think the non-conformance / error occurred? Use this area to explain your findings.

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>

<p>Assessment of the impact of root cause on patient results (Patient results may be affected before, during or after the PT event). Describe how the conclusion of impact was made and any corrective actions made to the patient result(s); "If No" - Explain why there was no patient impact?</p>
<p><b>IQC Values are within acceptable range.</b></p>
<p>Describe any previous proficiency testing issues with this test in the last sample:</p>

**Section 4: Department - Conclusion & Corrective / Preventive Action**

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>




<b>What corrective action have you carried out?</b>
Split Testing has been done , results are satisfactory.
<b>Preventive action put into place?</b>
NA
<b>Due date for closure of proposed corrective and preventive action:</b>
28-July-2024

**Person investigated**

**Department Manager**


**Lab Head**

**Signature with Date:**

Note: After signatures please hand over the form along with supporting documents to QA.

**Date when ANF received to QA along with supporting documents:** \_\_\_\_\_ **by:** \_\_\_\_\_

	<b>TATA 1mg Labs</b>	
	<b>TATA 1mg Technologies Private Limited</b>	
	<b>Form Name</b>	<b>Proficiency Testing - Action Needed Form</b>
	<b>Form No.</b>	<b>Gen / FR / 59</b>
	<b>Issue Date &amp; Version No.</b>	<b>01-Jul-2022 V2</b>