



Calibration Certificate

Calibration of the **B AUTO 200** (Fully Automatic chemistry Analyzer) has been carried out as per the recommendations of the principle company.

Customer Name : **STAR DIAGNOSTICS, ONGOLE, ANDHRA PRADESH**

Date of Installation : **20/02/2024 .**

Calibration Date : **16/07/2024.**

Next Calibration Due On: **15/01/2025.**

Tested By:

Checked By:

BEACON DIAGNOSTICS PVT. LTD.

FACTORY & REGISTERED OFFICE : 424, NEW GIDC, KABILPORE, NAVSARI - 396 424 (GUJ) INDIA
TOLL FREE : 1800-270-1027 | E-mail : customersupport@beaconindia.com | CIN No. : U24232GJ1990PTC013764

Performance Qualification

MODEL: B AUTO 200

SR.NO: AES2K0KH08055

Version: V4.8.3.7

Date: 16/07/2024

PQ is performed at

STAR DIAGNOSTICS, ONGOLE, ANDHRA PRADESH.

According to the assay results that are comprised in this document the B AUTO 200 systems can be approved for the routine operation.

PERFORMANCE CALIBRATION			
LOT No: 1228UE			
TEST	TARGET	RANGE	RESULT
GLU	270.00	230.00 ~ 310.00	294.83
AST	135.00	108.00~162.00	117.93
UREA	123.00	105.00~141.00	113.46
TBIL	4.95	3.91~5.99	4.32
DBIL	1.83	1.45~2.21	2.04
TP	4.78	3.82~5.74	4.04
ALB	3.25	2.76~3.74	3.60
CREA-E	5.06	4.06~6.06	4.40
TG	254.00	213.00~295.00	262.69
CHO	294.00	256.00~332.00	329.49

Note: Performance of the parameters are subject to traceability of Reagents, Manufacturer or Supplier.

Tested By:



Checked By:

BEACON DIAGNOSTICS PVT. LTD.

Operational Qualification

MODEL: B AUTO 200

Sr. No: AES2K0KH08055

Version: V4.8.3.7

Date: 16/07/2024

Operational Qualification Protocol run at STAR DIAGNOSTICS, ONGOLE, ANDHRA PRADESH

#	Description	Status
1	Blank/Optics	OK
2	Assays(Parameter) Calibration	OK
3	Assay (Parameter) Control – BEACON	OK
4	Assays(Parameter) Samples(Patient)	OK
5	Programming and operational training completed	OK
6	Instrument operating training complected	OK

I confirm that the above specified analyzer B AUTO 200 (Fully Automatic chemistry Analyzer) is installed properly and its performance is within the specification.

Tested By:


Checked By:

BEACON DIAGNOSTICS PVT. LTD.



BEACON GROUP



Installation Qualification

MODEL: BEACONIC B 200	Sr. No: AES2K0KH08055	Version: V4.8.3.7	Date: 20/02/2024
-----------------------	-----------------------	-------------------	------------------

Installation of the BEACONIC B 200 (Fully auto Chemistry Analyzer **STAR DIAGNOSTICS, ONGOLE, ANDHRA PRADESH**)
 The Installation Checklist is as given below:-

#	Position	Action	Status
1	Unpacking	Unpack the Complete System	OK
2	Inventory	Check Accessories With packing list	OK
3	Installation	Connection Of Pipeline Connection Of the Equipment Installation Of Software	OK
4	Service	Executing Air Purge Checking Syringe Sample and Reagent Probe Reaction Disk	OK
5	User	Start Work Procedure Check Reagent & DDW Cell blank Test Run QC Run Samples/Check Reproducibility Stop Work Procedure	OK

Tested By:

Checked By:

BEACON DIAGNOSTICS PVT. LTD.

FACTORY & REGISTERED OFFICE : 424, NEW GIDC, KABILPORE, NAVSARI - 396 424 (GUJ) INDIA
 TOLL FREE : 1800-270-1027 | E-mail : customersupport@beaconindia.com | CIN No. : U24232GJ1990PTC013764