



FACTORY & REGISTERED OFFICE :
424, NEW GIDC, KABILPORE, NAVSARI - 396 424. (INDIA)
PH. : +91-2637-265590, 265791 FAX : +91-2637-265105
E-mail : sales@beaconindia.com Website : www.beaconindia.com
CIN No. : U24232GJ1990PTC013764

Calibration Certificate

Calibration of the **B AUTO 400** (Fully Automatic chemistry Analyzer) has been carried out as per the recommendations of the principle company.

Customer Name : **ST. PIOUS HEALTHCARE PVT. LTD. BHAGALPUR, BIHAR**

Date of Installation : **28/11/2022.**

Calibration Date : **24/09/2024.**

Next Calibration Due On: **23/03/2025.**

KPN

Tested By:

Checked By:



FACTORY & REGISTERED OFFICE :
424, NEW GIDC, KABILPORE, NAVSARI - 396 424. (INDIA)
PH. : +91-2637-265590, 265791 FAX : +91-2637-265105
E-mail : sales@beaconindia.com Website : www.beaconindia.com
CIN No. : U24232GJ1990PTC013764

Performance Qualification

MODEL: B AUTO 400

SR.NO: AES2K1KH10040A

Version: V4.8.3.7

Date: 24-09-2024

PQ is performed at

ST. PIOUS HEALTHCARE PVT. LTD. BHAGALPUR, BIHAR.

According to the assay results that are comprised in this document the B AUTO 400 systems can be approved for the routine operation.

PERFORMANCE CALIBRATION			
LOT No: LP10-017			
TEST	TARGET	RANGE	RESULT
GLU	279.03	238.00 ~ 320.00	278.29
ALT	137.00	110.00 ~ 164.00	113.30
TBIL	6.06	4.00 ~ 6.12	5.05
ALB	2.45	2.51 ~ 3.39	3.14
UREA	120.00	102.00 ~ 138.00	121.77



FACTORY & REGISTERED OFFICE :
424, NEW GIDC, KABILPORE, NAVSARI - 396 424. (INDIA)
PH. : +91-2637-265590, 265791 FAX : +91-2637-265105
E-mail : sales@beaconindia.com Website : www.beaconindia.com
CIN No. : U24232GJ1990PTC013764

CHO	290.00	262.00 ~ 328.00	302.85
CREA-E	4.59	3.67 ~ 5.51	4.64
UA	9.31	8.09 ~ 10.53	9.42
TG	255.00	214.00 ~ 296.00	260.00

Note: Performance of the parameters are subject to traceability of Reagents, Manufacturer or Supplier.

Tested By:

Checked By:



FACTORY & REGISTERED OFFICE :
424, NEW GIDC, KABILPORE, NAVSARI - 396 424. (INDIA)
PH. : +91-2637-265590, 265791 FAX : +91-2637-265105
E-mail : sales@beaconindia.com Website : www.beaconindia.com
CIN No. : U24232GJ1990PTC013764

Operational Qualification

MODEL: B AUTO 400	Sr. No: AES2K1KH10040A	Version: V4.8.3.7	Date: 24-09-2024
-------------------	------------------------	-------------------	------------------

Operational Qualification Protocol run at ST. PIOUS HEALTHCARE PVT. LTD. BHAGALPUR,
BIHAR.

#	Description	Status
1	Blank/Optics	OK
2	Assays(Parameter) Calibration	OK
3	Assay (Parameter) Control – BEACON	OK
4	Assays(Parameter) Samples(Patient)	OK
5	Programming and operational training completed	OK
6	Instrument operating training completed	OK

- Error msg. list in case of hardware issue: Refer Operational Manual page no. 75

I confirm that the above specified analyzer B AUTO 400 (Fully Automatic chemistry Analyzer) is installed properly and its performance is within the specification.

Tested By:

Checked By:



FACTORY & REGISTERED OFFICE :
424, NEW GIDC, KABILPORE, NAVSARI - 396 424. (INDIA)
PH. : +91-2637-265590, 265791 FAX : +91-2637-265105
E-mail : sales@beaconindia.com Website : www.beaconindia.com
CIN No. : U24232GJ1990PTC013764

Installation Qualification

MODEL: BEACONIC B 400	Sr. No: AES2K1KH10040A	Version: V4.8.3.7	Date: 28-11-2022
-----------------------	------------------------	-------------------	------------------

Installation of the BEACONIC B 400 (Fully auto Chemistry Analyzer) at **ST. PIOUS HEALTHCARE PVT. LTD. BHAGALPUR, BIHAR.**

The Installation Checklist is as given below:-

#	Position	Action	Status
1	Unpacking	Unpack the Complete System	OK
2	Inventory	Check Accessories With packing list	OK
3	Installation	Connection Of Pipeline Connection Of the Equipment Installation Of Software	OK
4	Service	Executing Air Purge Checking Syringe Sample and Reagent Probe Reaction Disk	OK
5	User	Start Work Procedure Check Reagent & DDW Cell blank Test Run QC Run Samples/Check Reproducibility Stop Work Procedure	OK

Tested By:

Checked By: