CALIBRATION CERTIFICATE

CERTIFICATE NO: OCD/APL/0003 CALIBRATION DATE: 11/06/2024 CALIBRATED FOR: Lupin Diagnostics Limited CALIBRATION DUE: 10/12/2024

LOCATION : Chennai, Tamil Nadu. CALIBRATED AT: LABORATORY

SPECIFICATION OF INSTRUMENT:

INSTRUMENT	MODEL	SERIAL NO
VITROS CHEMISTRY	VITROS 250	25016302
System	J Number	J333425016302

Specification of Software:

SOFTWARE NAME	SOFTWARE VERSION	UPDATED ON
DOS Rmx os	V9.7	October 2019

The Reference of Calibration: The performance/Adjustments of various subsystems has been calibrated/tested by trained site engineer with pre-designed /calibrated tools provided for the particular subsystems by the company.

Preventive maintenance Procedure Checklists:

Subsystem Name

Adjusted/Verified

1. Sample Supply

* Hook Drive Belt

Software Adjustment:

* Tray Transport Hook Position

2. Slide Supply 1, 2 & slide transport

Mechanical Adjustment:

- * Dispense Blade Centering Adjustment
- * Dispense Blade Drive Belt Tension
- * Plunger Centering Alignment

Software Adjustment:

* Dispense To Metering Position

3. Sample Metering

Mechanical Adjustment:

* Truck Sensor

√

Ortho Clinical Diagnostics India Private Limited
Regd office: 403,Leela Business Park,Andheri Kurla Road ,Andheri(E),Mumbai 400059.
T:+91 22 6787 9300 F:+91 22 6787 9333 CIN: U51397MH2015FTC262650

Email ID: connectocd@orthoclinicaldiagnostics.com Website: www.orthoclinicaldiagnostics.com

Ortho Clinical Diagnostics

- * Truck Guide Button
- * Metering Truck
- 4. Software Adjustment:
 - * Tip Pickup X, Y and Height
 - * Sample Tray at Sample Position –Auto Sample Tray
 - * Proboscis Centering and Sample Cup Depth
 - * Mixing Cup X, Y and its Depth
 - * Diluents Tray
 - * Rate / Colorimetric Metering Position
 - * Potentiometric Metering Position
 - * Proboscis Compression
 - * Wetness Detector
- 5. Slide Incubator

Software Adjustment:

- * Incubator Disk Sizing
- * Incubator Disk Stopping
- * PM Slide Insertion Depth
- * Slide Read Position (RT/CM Depth & Read Sync)
- * Potentiometric Eject Blade Depth
- 6. Reflectometer
 - * IRIS Lamp Output Voltage
 - * Correction Factor
- 7. Master Computer
 - * Touch Screen Calibration
 - * System Full Backup

CALLIBRATED BY



Naveen K – Service Engineer (CTS)

- 1. This certificate refers only to the particular item submitted for calibration.
- 2. The calibration result reported in the certificate is valid at the time of and under the stated condition of the measurement.

Website: www.orthoclinicaldiagnostics.com