



SRI SELLAM
Clinical Lab

113/4, Sulthanpet
Near Panchamuga Vinayagar Kovil
P. Velur-638182, Namakkal Dt.
Off : 04268-222750
Off : 7373734265
E-mail: srisellamlab@gmail.com

R.Muthusamy B.Sc., Bio Chemist

Mobile : 94423 34265, 73730 34265

REPORT
TO WHOMSOEVER IT MAY CONCERN

Dear Madam,

We will apply the EQAS for ISHTM-AIIMS. As we are requested to our ILC lab to provide us the report, they are reverted us it will be given after a week because of COVID19 Situation, once as we received the Document sure we will share the report of comparison through mail. Kindly accept our request and make us favour.

Thank you

Yours Sincerely

Muthusamy R
Proprietor
Sri Sellam Clinical Lab

SRI SELLAM CLINICAL LAB
113/4, Sulthanpet
Near Panchamuga Vinayagar Kovil
P. VELUR-638182, Namakkal Dt.
Cell: 94423 34265



ISHTM-AIIMS EXTERNAL QUALITY ASSURANCE PROGRAMME

NABL accredited programme as per ISO/IEC 17043:2010 standard

Organized by

Department of Hematology, AIIMS, New Delhi-110029

E-Mail ID: accuracy2000@outlook.com , Contact No. 9013085730



New Registration Confirmation

To ,SRI SELLAM CLINICAL LAB 113/4, SULTHANPET, NEAR PANCHAMUGA VINAYAGAR KOVIL, PARAMATHI VELUR-638182

Dated: Nov 19,2020

Payment Detail:

Amount	Transaction No	Transaction Date:
2009.44	TXN1605765314	2020-11-19

** This is electronic copy for the confirmation of payment submission to ISHTM-AIIMS-EQAP and does not need any signature.

** Average time for allotment of the EQAP Code No. and first sample receipt takes an average time of 2-3 months, please wait.

-----End Of Receipt-----

With Regards.

Dr. Seema Tyagi (Prof, Hematology)

ISHTM-AIIMS-External Quality Assurance Programme

(An ISO 17043:2010 certified programme)

Room No. 209, Department of Hematology

2nd Floor New Pvt Ward

All India Institute of Medical Sciences

New Delhi-110029

Website: www.ishtmaiimseqap.com

This is an automatically generated e-mail, please do not reply to it. For any further enquiry e-mail us only on accuracy2000@gmail.com (ISHTM-AIIMS-EQAP)



© ISHTM-AIIMS-EQAP