

Original Buyer Copy

TAX INVOICE

[See Rule 46 of the Central Goods and Services Tax Rules, 2017]

TULIP DIAGNOSTICS (P) LTDRegd Off. "Gitanjali"-Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex
PO. Goa 403 202 INDIA.

IRN: 67a04c24963ef45511d729565c37712e2057542d66280631c5ad97db1bb79656

INV NO: TAI1/0256/2021

29-Jan-2021

Supplier Details

Name TULIP DIAGNOSTICS (P) LTD
 Division Name
 Regd Office "Gitanjali"-Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O. Goa 403 202 INDIA
 GSTIN Unique ID 30AAACT1411B1ZH
 PAN NO AAACT1411B
 State Goa
 Works Address
 Warehouse Address
 Ref Date
 Mfg DL Nos N/A
 Wholesale DL Nos NA
 Transport

Customer Details / Bill To

Code C23B0009
 Customer Name BANSAL HOSPITAL
 Customer Billing C-SECTOR, SHAHPURA, BHOPAL (M.P)
 Address
 Billing District/ State Madhya Pradesh
 PAN NO
 Billing GSTIN 23AABCA2560P1ZP
 Billing Mob No Shipping Mob No
 Billing State Madhya Pradesh
 Service To
 Consignee BANSAL HOSPITAL
 Name
 Shipping Address C-SECTOR, SHAHPURA, BHOPAL (M.P)
 Shipping Address Madhya Pradesh 23AABCA2560P1ZP
 District/ State & Code

Sr.No	Product Description	UOM	Qty	Rate	Taxable Value	SGST		CGST		IGST		Total Amount
						Rate %	Amnt ₹	Rate %	Amnt ₹	Rate %	Amnt ₹	
			20-Jan-21	TO	19-Jan-22							
1	ANNUAL MAINTENANCE CONTRACT SAC 998719 For: Lisawash - 4000Sr.No. 591602027-E3A AMC Period From 20-Jan-21 TO 19-Jan-22	NOS	1	11000.0000	11000.00	0.00	0.00	18%	1980.00	0	12980.00	
2	ANNUAL MAINTENANCE CONTRACT SAC 998719 For: Lisaquant - TSSr.No. 451527026E AMC Period From 20-Jan-21 TO 19-Jan-22	NOS	1	11000.0000	11000.00	0.00	0.00	18%	1980.00	0	12980.00	

Remarks

Note

Declaration: We declare that the supplier price shown above is sole consideration and certified that particulars given above are true and correct and the amount indicated represents the price actually charged and that there is no flow of additional consideration directly or indirectly from the buyer.

Amount in Words: Rupees Twenty Five Thousand Nine Hundred and Sixty Only

TAXABLE AMOUNT ₹	22000.00
SGST ₹	0.00
CGST ₹	0.00
IGST ₹	3960.00
NET AMOUNT ₹	25960.00

For, TULIP DIAGNOSTICS (P) LTD



Installation cum Warranty Certificate

Name & Address of Customer:

BANSAL HOSPITAL
C - Sector, Shahpura
Bhopal (M P)

Instrument Model: LISAWASH - 4000

Instrument Sr. No.: 591602027-E3A

Installation Date: 19/05/2014

Invoice No. TR/019/15

Invoice Date: 24/04/2014

AMC Due Date 18/05/2016

Contact Person: Mr. Vikash Gupta.
Phone No.: 9893845512

Power supply details:

Instrument connected to: Online UPS Yes/No Offline UPS Yes/No

Supply input to Machine from UPS: LN: NE: LE:

Supply input to UPS from Mains: LN: NE: LE:

Accessories Provided: main system, 8-pin manifold, Bottles, (a) waste bottle, (b) wash bottle (c) water bottle, power cable, user manual, fuses, Dust cover, cleaning wipe, operation procedure card.

Persons Trained:

Mr. Vikash Gupta

Demonstrated:

yes.

Remarks: If any

Note: The above machine is satisfactorily installed at the site and staff adequately trained
Received all accessories as per packing list



Name of Department Head:

Signature & Stamp

Date:

Company's Sales & Service Engineer:

Signature

Date:

[Signature]
19/05/14



[Signature] Vikash Gupta
Tulip Diagnostics P Ltd.

Gatani, Tulip Block, Dr. Antonio De Rego Bagh, Alto Santacruz, Bambolim P.O. Goa. Pin: 403 202.
Tel: 0832-2455400, Fax: 0832-466444, E-mail: sales@tulipgroup.com, Website: www.tulipgroup.com

Lic.No.-28C/3/2014
Bansal Hospital, Bhopal

TULIP DIAGNOSTICS (P). LTD. (INSTRUMENT DIVISION)

BAMBOLIM, ST-INEZ & VERNA - GOA

Ref. SOP No: SOP/ID/013
CC No: CC-18/ID/001

CALIBRATION CERTIFICATE

Model: Lisawash 4000

Sr. No.: 591602027 E3A

Calibration Date: 25/01/2021

Calibration Due Date: 24/01/2022

Check Base plate & Top Cover Visually

Tick yes if OK, or else tick no.

- | | | |
|---|---|-----------------------------|
| 1) Prime: - Level in the rinsing strip | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| 2) Rinse: - Level in the rinsing strip | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| 3) Plate wash: - Residual Volume (<3 µl) | | |
| Plate no. 1 400µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| Plate no.2 350µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| Plate no.3 300µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| 4) Strip wash: - Residual Volume (<3 µl) | | |
| Plate no. 1 400µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| Plate no.2 350µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| Plate no.3 300µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| 5) Dispensing: - Dispensing is without splashing or not | | |
| Without Rinsing | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| With Rinsing | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| Dispensing Volume: | | |
| Plate no. 1 400µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| Plate no.2 350µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| Plate no.3 300µl | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| 6) Waste bottle: - Check for Waste Bottle sensor | | |
| Message on Display & Beep when Bottle is Full | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |
| 7) Wash1 bottle: - Check for Wash1 Bottle sensor | | |
| Message on Display & Beep when Bottle is empty | YES <input checked="" type="checkbox"/> | NO <input type="checkbox"/> |

FM/ID/029-00

Page 1 of 2

Corporate Office: Gitarjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim Complex P.O., Goa - 403 202, INDIA.
Tel.: (0832) 2458546-51, Fax: (0832) 2458544, E-mail: sales@tulipgroup.com, Website: www.tulipgroup.com

Innovative Excellence

TULIP DIAGNOSTICS (P). LTD. (INSTRUMENT DIVISION)

BAMBOLIM, ST-INEZ & VERNA - GOA

Ref. SOP No: SOP/ID/013
CC No: CC18/ID/001

- 8) Wash2 bottle: - Check for Wash2 Bottle sensor
Message on Display & Beep when Bottle is empty YES NO
- 9) Wash3 bottle: - Check for Wash3 Bottle sensor
Message on Display & Beep when Bottle is empty YES NO
- 10) Water bottle: - Check for Water Bottle sensor
Message on Display & Beep when Bottle is empty YES NO
- 11) Sensor cable: - Check for Sensor cable
Beep when sensor cable is not connected YES NO
- 12) Incubator: - Check for temperature
Ambient, 25 °, 30 °, 37° C ± 1° C YES NO

Measured Temperature	25°C, 30°C, 37°C
Displayed Temperature	25°C, 30°C, 37°C

13) Tools used

Description	Sr. No/ Model	Traceability
Digital Multimeter	MASTECH/MS8264	QATECH/GOA/01/20/IN121000002

14) Keyboard: -

Check all keys on keyboard by pressing
Them individually, beep should be heard when pressed. YES NO

15) Display

Check all the characters are visible properly. YES NO

16) Horizontal Movement

There should not be any sound or vibration YES NO

17) Vertical Movement

There should not be any sound or vibration. YES NO

18) Centering check

Check by reversing the Plate YES NO

OTHER INSTRUCTIONS (IF ANY, PLEASE SPECIFY)

CALIBRATED BY:

ARM



APPROVED BY:

[Signature]
Eloyd Dsouza

FM/ID/029-00

Page 2 of 2

Installation cum Warranty Certificate

Name & Address of Customer

BANSAL HOSPITAL
C-Sector Shahpura, Bhopal (M.P)

Instrument Model: LISAQUANT - TS

Instrument Sr. No.: 451527026E

Installation Date: 19/05/2014

Invoice No. TR/1189/14

Invoice Date: 25/03/2014

Contact Person: Mr. Vikash Gupta

Phone No: 9593845512

AMC Due Date: 18/05/2016

Power supply details

Instrument connected to Online UPS Yes/No Offline UPS Yes/No

Supply input to Machine from UPS LN NE: LE:

Supply input to UPS from Mains LN NE: LE:

Accessories Provided: Operator Manual, Mouse, Jugs, Dust cover, power cord, Touch pen & pen Bracket.

Persons Trained: Mr. Vikash Gupta

Demonstrated: Yes.

Remarks if any

Note: The above machine is satisfactorily installed at the site and staff adequately trained
Received all accessories as per packing list



Name of Department Head:

Signature & Stamp

Date:

Company's Sales & Service Engineer:

Signature

Date:



Vikash Gupta

19/05/14

Tulip Diagnostics P Ltd.

Gitanjali, Tulip Block, Dr. Antonio Do Rego Bagh, Alto Santacruz, Bambolim P.O. Goa. Pin: 403 202.

Call: 0832-2458546/51 Fax: 0832-2458544, E-mail: sales@tulipgroup.com, Website: www.tulipgroup.com

BLOOD BANK & COMPONENT CENTER
Lic.No.-28C/3/2014

Bansal Hospital, Bhopal

TULIP DIAGNOSTICS (P). LTD. (INSTRUMENT DIVISION)

BAMBOLIM, ST-INEZ & VERNA - GOA

Ref. SOP No: SOP/ID/013

CC No: CC18/ID/001

CALIBRATION CERTIFICATE

Model: Lisaquant - TS

Sr. No.: 451527026E

Calibration Date: 25/01/2021

Calibration Due Date: 24/01/2022

Check Base plate & Top Cover Visually

1) A/D VALUES

- WITHOUT MICROWELLS (AIR) (Range 6000-65000)

Tick YES if OK, or else tick NO.

FILTER	405	450	492	630
A	35737	36312	35810	36410
B	35182	35918	35318	36513
C	35768	35810	35751	36151
D	35124	36185	35789	36312
E	35218	36315	36011	36189
F	35310	36213	35721	36700
G	35217	36189	35719	37000
H	35157	35713	36319	36011

With in range

YES

NO

2) TOUCH PANEL

Check the all Menu by clicking them individually with help of Touch pen, beep should be heard when clicked.

YES

NO

3) DISPLAY

Check all the Menu & characters on display are visible properly.

YES

NO

4) EXTERNAL PRINTER CONNECTION

Check printing.
All characters should print properly.

YES

NO

5) Horizontal Movement

There should not be any sound or vibration.

YES

NO

6) Centering checked

By reversing the filmstrip check readings are Repeating or not.

YES

NO

FM/ID/029-00

Page 1 of 2

LIP DIAGNOSTICS (P). LTD. (INSTRUMENT DIVISION)

BAMBOLIM, ST-INEZ & VERNA - GOA

Ref. SOP No: SOP/ID/013
CC No: CC18/ID/001

7) MONOCHROMATIC DISTILLED WATER READINGS

YES

NO

Check water readings in Absorbance mode. It has to be Below 0.150 in new plate with 150µl water.

Reading of distilled water 450nm filter.

A1 well	B1 well	C1 well	D1 well	E1 well	F1 well	G1 well	H1 well

8) BIOCHROMATIC DISTILLED WATER READINGS

YES

NO

Check water readings in Absorbance mode. It has to be above 0.000 & below 0.050 in new plate with 150µl water.

Reading of distilled water 450nm and 630nm filter.

A1 well	B1 well	C1 well	D1 well	E1 well	F1 well	G1 well	H1 well

PASSED

(Tick yes if OK or else tick no, and give reason)

YES

No

CALIBRATED BY:

Handwritten signature



APPROVED BY:

Handwritten signature
[Floyd Osaura]

FM/ID/029-00

Page 2 of 2