



CHRISTIAN MEDICAL COLLEGE

DEPARTMENT OF CLINICAL BIOCHEMISTRY

CMC EXTERNAL QUALITY ASSURANCE SCHEME

MONTHLY SUMMARY REPORT - JULY 2022



PC-1024

Lab Name MDC X RAYS LABORATORY & ECG

Lab No

15334

Constituent Group Chemistry I

Date of Result Entered :

11/07/2022

PT item Lyophilized Serum

Date of Report Published :

05/08/2022

Sl.No	Analyte	Method / Principle Name	Analyzer Name	No of Participants	DV	Participants		Your Value	SDI	U
						CV	SD			
1	GLUCOSE	Hexokinase	Roche C111, C311, Integra 400 plus, C501, C502, C702	293	308.62	3.96	12.21	298.7 mg/dL	-0.81	1.43
2	UREA	Urease UV / GLDH	Roche C111, C311, Integra 400 plus, C501, C502, C702	300	129.16	5.58	7.20	109.2 mg/dL	-2.77	0.83
3	CREATININE	Jaffes Kinetic - Alkaline picrate	Roche C111, C311, Integra 400 plus, C501, C502, C702	256	4.48	5.49	0.25	4.15 mg/dL	-1.34	0.03
4	T.BILIRUBIN	Diazonium salt (Colorimetric) / Jendrassik	Roche C111, C311, Integra 400 plus, C501, C502, C702	292	0.95	7.14	0.07	0.9 mg/dL	-0.74	0.01
5	T-PROTEIN	Biuret - Colorimetric	Roche C111, C311, Integra 400 plus, C501, C502, C702	304	5.37	4.25	0.23	4.8 g/dL	-2.50	0.03
6	ALBUMIN	BCG - colorimetric	Roche C111, C311, Integra 400 plus, C501, C502, C702	292	3.50	4.69	0.16	3.35 g/dL	-0.91	0.02
7	CALCIUM	5-Nitro-5-methyl-BAPTA	Roche C111, C311, Integra 400 plus, C501, C502, C702	236	9.99	4.08	0.41	9.2 mg/dL	-1.94	0.05
8	PHOSPHORUS	Molybdate UV / Phosphomolybdate complex	Roche C111, C311, Integra 400 plus, C501, C502, C702	235	6.85	4.33	0.30	6.3 mg/dL	-1.85	0.04
9	URIC ACID	Enzymatic / Uricase Colorimetric	Roche C111, C311, Integra 400 plus, C501, C502, C702	292	3.04	4.61	0.14	3 mg/dL	-0.29	0.02
10	CHOLESTEROL	CHOD-PAP	Roche C111, C311, Integra 400 plus, C501, C502, C702	308	118.77	4.35	5.17	106.6 mg/dL	-2.35	0.59
11	TRIGLYCERIDE	GPO-PAP / Enzymatic Colorimetric / End Point	Roche C111, C311, Integra 400 plus, C501, C502, C702	277	273.68	4.16	11.38	249.2 mg/dL	-2.15	1.37
12	HDL	Direct method / Enzymatic colorimetric	Roche C111, C311, Integra 400 plus, C501, C502, C702	273	23.20	4.87	1.13	21.9 mg/dL	-1.15	0.14
13	SODIUM	ISE - Direct	Any Analyser (Automation / Semi Automation)	1047	139.31	3.44	4.80	130.1 mmol/L	-1.92	0.30
14	POTASSIUM	ISE - Direct	Any Analyser (Automation / Semi Automation)	1099	5.06	6.41	0.32	4.59 mmol/L	-1.45	0.02
15	CHLORIDE	ISE - Direct	Any Analyser (Automation / Semi Automation)	878	100.12	4.63	4.64	99.9 mmol/L	-0.05	0.31
16	AST	UV kinetic(with & without PLP (P-5-P))	Roche C111, C311, Integra 400 plus, C501, C502, C702	266	124.86	6.12	7.65	119.6 U/L	-0.69	0.94
17	ALT	UV kinetic(with & without PLP (P-5-P))	Roche C111, C311, Integra 400 plus, C501, C502, C702	266	55.82	9.22	5.15	54.5 U/L	-0.26	0.63
18	ALP	PNP AMP kinetic	Roche C111, C311, Integra 400 plus, C501, C502, C702	281	211.42	7.02	14.85	203.6 U/L	-0.53	1.77

SDI Range	Interpretation
Within -1.00 to +1.00	Excellent.
Within ± 1.01 to ± 2.00	Good.
Within ± 2.01 to ± 2.99	Accept with caution. Warning Signal.

LAB ADDRESS :
 MDC X RAYS LABORATORY & ECG
 HEALTH CLUB BUILDING MANANTHAVADY
 WAYANAD
 KERALA670645

Coordinator Contact Details:
 Email: clinqc@cmcvellore.ac.in
 Contact Number: 0416-2283102

Pamela Christudoss
 Dr. Pamela Christudoss
 CMC EQAS Coordinator
 Christian Medical College, Vellore

Homogeneity and Stability of the sample is passed.
 Data in CMC EQAS reports is confidential
 CMC EQAS does not sub contract any components
 ***** End of Report *****

S.L NO.	Assay	Corrective Action Taken	QC value Observed	Remarks	Sign.
1.	Urea	Recalibrated	Within range.	QC Passed.	<i>B</i>
2.	Total Protein	Recalibrated	Within range.	QC Passed	<i>B</i>