



Beckman Coulter India Pvt.

Regional Office

The DLF Galleria bearing unit no. 1201,
DLF Galleria, Plot No. BG-08,
Premises No. 02-0124, Action Area IB,
New Town, Kolkata – 700156.

Tel: +91 33-40050630/2

Fax: +91 33-40050631

CALIBRATION CERTIFICATE

We here by certify that the fully automated chemistry analyzer Beckman Coulter AU480 has been installed by us and Calibrated. It has passed all the following quality checks.

Instrument serial number : 2018041921

Installed at : Sri Venkateswara Institute Of Medical Sciences
Tirupati

Last Calibration Done on : 20.04.2022

Next Calibration Due on : 19.04.2023

Quality Checks - Passed

- . Reagent and sample aspiration system.
- . Reagent Tray & STAT Table Refrigeration.
- . Cuvette Wheel Incubation.
- . Lamp & Optical Unit.
- . Water supply and Tubing system.
- . Calibration of Probes, Rack, Track & Carousel Movements.
- . Reagent and Sample Position Identification.
- . Water test and QC checks.
- . Calibration of all test parameters.
- . QC run for all test parameters.

We here by confirm that AU480 has passed all the above quality checks and is working in good condition.

Beckman Coulter India Pvt. Ltd

Registered Office: Unit Nos. TF-B-07(A) to 15, Third Floor, B Wing, Art Guild House, Phoenix Market City, LBS Road,
Kurla (West), Mumbai – 400 070. Tel: +91 22 3080 5000 Fax: +91 22 3080 5060
CIN: U33119MH2005PTC157177 www.beckmancoulter.com

Confidential - Company Proprietary



Beckman Coulter India Pvt.

Regional Office
The DLF Galleria bearing unit no. 1201,
DLF Galleria, Plot No. BG-08,
Premises No. 02-0124, Action Area IB,
New Town, Kolkata – 700156.
Tel: +91 33-40050630/2
Fax: +91 33-40050631

CALIBRATION CERTIFICATE

Instrument serial number : 2018041921
Installed at : Sri Venkateswara Institute Of Medical Sciences
Tirupati

Last Calibration Done on : 20.04.2022
Next Calibration Due on : 19.04.2023

Photometric check

Requirement : All O.Ds should be less than 1.5

Wavelength	Optical Density	Result
340	1.2385	OK
380	0.8929	OK
410	0.8776	OK
450	0.8664	OK
480	0.7470	OK
520	0.7611	OK
540	0.8359	OK
570	0.9501	OK
600	1.0537	OK
660	0.8075	OK
700	0.8017	OK
750	0.6111	OK
800	0.6009	OK

Beckman Coulter India Pvt. Ltd

Registered Office: Unit Nos. TF-B-07(A) to 15, Third Floor, B Wing, Art Guild House, Phoenix Market City, LBS Road,
Kurla (West), Mumbai – 400 070. Tel: +91 22 3080 5000 Fax: +91 22 3080 5060
CIN: U33119MH2005PTC157177 www.beckmancoulter.com

Confidential - Company Proprietary



Beckman Coulter India Pvt.

Regional Office

The DLF Galleria bearing unit no. 1201,
DLF Galleria, Plot No. BG-08,
Premises No. 02-0124, Action Area IB,
New Town, Kolkata – 700156.

Tel: +91 33-40050630/2

Fax: +91 33-40050631

Incubation Temperture

Range: 37C +0.5

On display the temperature is 37.0 stable

S.NO	TEST PARAMETERS	MEASURED/CALIBRATED	REMARKS
01	INPUT VOLTAGE	230 V	--ok--
02	POWER	4 KVA	--ok--
03	INPUT FREQUENCY	50 Hz	--ok--
04	OPERATING TEMP.	20-25 °C	--ok--

Tool Used: Digital Multimeter

Make: Fluke

Model: 17B+

Sl.no: 39591421WS

We hereby confirm that Beckman Coulter AU480 has been successfully calibrated and passed all the above quality checks and is working in good condition.

For Beckman Coulter India Pvt Ltd

**T V SAPTHAGIRIDHAR
CUSTOMER SUPPORT ENGINEER**

Beckman Coulter India Pvt. Ltd

Registered Office: Unit Nos. TF-B-07(A) to 15, Third Floor, B Wing, Art Guild House, Phoenix Market City, LBS Road, Kurla (West), Mumbai – 400 070. Tel: +91 22 3080 5000 Fax: +91 22 3080 5060
CIN: U33119MH2005PTC157177 www.beckmancoulter.com

Confidential - Company Proprietary