

2 & 3, Devi Illam, Priya Nagar, Telungupalayam Pirivu, Coimbatore - 641 039, TamilNadu

T. ARUN PRASAD, M. Tech., MBA., (Ph.D)

Managing Director
Consultant Biomedical Engineer

GSTIN: 33APLPA1828A1ZM

## **Calibration Protocol & Report for Medical Equipment**

Equipment: Micro pipette

Protect Class/Type:

B

Make

: micro lux

SNo

NM/MP/001

Mode!

: 100-1000ul

Location

**KOTHAGONDAPALLI - GPHC** 

Measurement Parameters			Status	Remarks
Electrical test	(a). oral r	(a). oral method digital apparatus		well
II Physical Condition Test		(a). Cabinet & Visible damages, if any (b). on/off switch, digital display,glass		Checked ok
		tube visual checkups (c). Risk of mechanical damages, if any		-
Performance		(a). Glass visual clarity		Spares
Test	(b). oral sensor edge		ok	replacement/
	(c).Abnormalities		No	Checked &
	(d). Measurement checkup		Ok	measured ok
	(e). on/off functioning		Well	
	(f). devia	(f). deviations, if any		
IV Calibration Test	(a). digital zero calibration tested with		zero level	ensured zero
	i. zero degree (with freezer)		confirmed/	level. fine
	ii. normal body temp		Checked ok	
	(w:	(warmer/incubator)		
	iii. 100	iii. 1000micro ml (water)		
	(b). Meas	(b). Measurement correlation test		
Formulated By: Einstein Research Issued By: Laboratory		User Department:		
sed By :	<b>A</b>	SILICON TECHNOLOGIES	. 1	
ved On : 12.09.23		COIMBATORE - 641 039	<b>7</b> !	
eview : 12.09.24		Cell: 98652 85754, 81246 75020.	1	
	Physical Condition Test  Performance Test  Calibration Test  Idated By: Einstein Story Sed By: Sch team head Syd On: 12.09.23	Physical Condition Test  Physical Condition Test  Performance Test  Performance (a). Glass (b). oral s (c). Abnor (d). Meas (e). on/of (f). devia  i. zer ii. no (w. iii. 100 (b). Meas  lated By: Einstein Research story sed By: ch team head ved On: 12.09.23	Electrical test  (a). oral method digital apparatus  Physical Condition Test  (b). on/off switch, digital display, glass tube visual checkups (c). Risk of mechanical damages, if any  Performance Test  (a). Glass visual clarity (b). oral sensor edge (c). Abnormalities (d). Measurement checkup (e). on/off functioning (f). deviations, if any  Calibration Test  (a). digital zero calibration tested with i. zero degree (with freezer) ii. normal body temp (warmer/incubator) iii. 1000micro ml (water) (b). Measurement correlation test  Issued By: SILICON TECHNOLOGIES A Biomedical Instrumentaria Companiculation on paniculation companiculation companiculati	Electrical test  (a). oral method digital apparatus  (b). on/off switch, digital display, glass (c). Risk of mechanical damages, if any (d). oral sensor edge (c). Abnormalities (d). Measurement checkup (e). on/off functioning (f). deviations, if any  Calibration Test  (a). digital zero calibration tested with i. zero degree (with freezer) ii. normal body temp (warmer/incubator) iii. 1000micro ml (water) (b). Measurement correlation test  Issued By: Checked ok  Well  SILICON TECHNOLOGIES A Biomedical fist unless against a company COMMEDICAE - 641 039  Ok  No  Ved On: 12.09.23  Ok  No  No  No  No  SILICON TECHNOLOGIES A Biomedical fist unless against a company COMMEDICAE - 641 039

Type your text

