Ortho Clinical Diagnostics

403, Leela Business Park, Andheri Kurla Road, Andheri East, Mumbai – 400059

T:+91 22 6787 9300 F:+91 22 6787 9333

CALIBRATION CERTIFICATE

CERTIFICATE NO:OCD/ESICHOSPITAL/003

CALIBRATED FOR :ESICHOSPITAL, ERRAGADDA.

LOCATION

:MICROBIOLOGY

CALIBRATION DATE: 23/12/2022

CALIBRATION DUE: 22/12/2023

CALIBRATED

: AT SITE

SPECIFICATION OF INSTRUMENT:

INSTRUMENT	MODEL	SERIAL NO	
Vitros	VITROS ECI	3000-6046	
IMMUNODIAGNOSTICS SYSTEM	J NUMBER	J3000-6046	

Specification of Software:

SOFTWARE NAME	SOFTWARE VERSION
QNX os	V3.9.2

The Reference of Calibration: The performance of various subsystems has been calibrated/tested by trained site engineer with pre-designed calibrated tools provided for the particular subsystems by the company.

Subsystem 1.Sample Sup	Name Verified ply Inspection	
2.Sample Met		
	* SampleTrayTipPickup & SamplingCenter	
	* Auxiliary Tip Supply * Sample Dispense	
	* Outer Ring Mapping	
2 Descent Cur	* Sample Metering Leak Test	
3.Reagent Sup		
e 2 -	* Carousel to Well Shuttle	
	* Well Shuttle – Positioning	
	* Pack Opener (Home/Open/Close)	
	* Outer Reagent – Alignment	
	* Inner Reagent – Alignment	
	* Well Top Alignment	\neg
4.Reagent Me		
	* Wash Station	
w	* Dispense	
5.Processing (Cente	
	* Inner Lift Pin (Home/Shuttle/Read)	
	* Inner/Outer rings Home Position	
	* Incubator Shuttle	

* Incubator Temperature calibration

6. Well Wash Assembly

Ortho Clinical Diagnostics

403, Leela Business Park, Andheri Kurla Road, Andheri East, Mumbai – 400059 T:+91 22 6787 9300 F:+91 22 6787 9333

* We	11 W	Vash	Probe	Horizo	ntal	Alignme	ent
------	------	-------------	-------	--------	------	---------	-----

* Inner Ring Mapping

* Station Verification

7. Signal Reagent Assembly

* Signal Reagent Probe - Horizontal Alignment

* SR Dispense

* SR Station Verification

8.Luminometer

* Full Calibration (optional)

* IRS Calibration



Pressure Regulator Calibration

* Supply Level Calibration







The results of comparison are as follows:

SUBSYSTEM NAME	CALIBRATION ACCURACY		CORRECTION	
	RANGE	OBTAINED		
SAMPLE METERING	AUTOMATIC			
REAGENT METERING	4.90ML-5.10ML	5.0	ML	
DISPENSE CALIBRATION				
REAGENT METERING	2.05ML-2.55ML	2.3	ML	
OUTER PROBE WASH		x " 12 =		
WELL WASH	1.8ML-2.4ML	2.0	ML	
ASPIRATION				
WELL WASH DISPENSE	6.2ML-6.4ML	6.3	ML	
SIGNAL REAGENT	9.99ML-10.01ML		ML	
DISPENSE PUMP A		10.0		
SIGNAL REAGENT	9.99ML-10.01ML	10.0	ML	
DISPENSE PUMP B		7		

CALIBRATED BY



Anuganesh Kasthuri Sr Service Engineer Ortho care srvice Ortho Clinical Diagnostics

NOTE:

- 1. This certificate refers only to the particular item submitted for calibration.
- 2. The calibration result reported in the certificate is valid at the time of and under the stated condition of the measurement.