

Product Name	Auto analyzer	Model SN Standard and Demand	CSENSE 200 ^{CD} 0532405005 IEC		
fested Date	02.08.2023				
No.			Obtained	Conclusion	
ī	Basic Function	Effective	Working Satisfactory	Ø	QC pass
2	Calibration and QC programming	Effective	Done	8	QC pass
3	Sample entry	Effective	Done	8	QC pass
4	All wavelength readings	20k to 50k	Working Satisfactory	V	QC pass
5	Reaction well temperature	37°C ± 0.3	37.1°C	∞	QC pass
6	Reagent cooling	15°C lesser than RT	Checked	8	QC pass
7	Functionality check of mixer motor.	Judge whether the mixer is in the centre of the test cuvette/cleaning pool and whether you can hear the sound of rubbing up against the wall at the time of rotation	Checked	8	QC pass
8	Mechanical co-ordinate values	Instrument specific	Checked	V	QC pass
,	level sensor	Effective	Checked	8	QC pass
10	Collision sensor	Effective	Checked	8	QC pass
11	Input power source	210 to 240 V,50Hz	232 V, 50HZ	V.	QC pass
Testing Result :	Passed alibrated as per the factory standard	QC PASS and the perfromed check values are satisfac	tory.	18/ 1	Stomedicals 3
Next calib	oration on: 01.08.2024			SSO CE	ned :

Medsource Ozone Biomedicals Pvt. Ltd.