




Date: 14-1-2024

TO WHOSOEVER IT MAY CONCERN

The medica electrolyte analyser NA/K/CL Sr.no 63557cnkc placed at Vinamra Swaraj Hospital Vashi is under AMC Contract of Transasia Bimedicals LTD.

AMC end date : 13-1-2025

from 14-1-2024 to 13-1-2025 is under annual maintaince contract.



Thanks & Regards,

Suraj Chandak

ASM

9136950458



UNMATCHED SERVICE
SINCE 1979....

02-04-2017

TO WHOMSOEVER IT MAY CONCERN

THIS IS TO STATE THAT EASYLITE PLUS DOES NOT REQUIRE
EQUIPMENT CALIBRATION.



Calibration Certificate



UNMATCHED SERVICE
SINCE 1979

The EasyLyte displays analysis results automatically. It is a Self Calibrated analyzer. It is programmed to Calibrate every four hours.

Calibration is automatic, but can be performed on demand. A unique solutions pack contains standard solutions to calibrate the analyzer. The disposable solutions pack also contains a waste receptacle for spent samples and solutions, eliminating the need to handle biological waste materials.

Calibration values are the electrode response values (slopes) during the last calibration. Acceptable ranges are as follows:

Na+: 50-64

K+: 50-64

Cl-: 40-64

Li+: 50-66

Ca++: 20-34

pH: 45-64

MEDICA EasyLyte